

ColorFlex® EZ

Supplemental Manual for EasyMatch® QC



Hunter Associates Laboratory
11491 Sunset Hills Road
Reston, Virginia 20190 USA
www.hunterlab.com

A60-1017-662
Manual Version 1.1

Copyrights and Trademarks

This documentation contains proprietary information of Hunter Associates Laboratory, Inc. Its reproduction, in whole or in part, without express written consent of Hunter Associates Laboratory, Inc. is prohibited.

EasyMatch QC and ColorFlex are registered trademarks for Hunter Associates Laboratory, Inc.

Windows is a registered trademark of Microsoft Corporation in the United States and other countries.

Duraflect, Spectrafect, and Spectralon are trademarks of Labsphere, Inc.

Teflon is a registered trademark of Dupont.



Caution: If the equipment is used in a manner not specified by the HunterLab, the overall safety may be impaired. The instrument is for indoor use only and not suitable for a wet location.



Caution: There is a potential of a UV Light hazard in using this instrument. Please avoid looking directly at the light.

Contents

COLORFLEX EZ FEATURES	5
ColorFlex EZ Accessories	6
COLORFLEX EZ INSTALLATION	7
Install EasyMatch QC Software	7
Activate the SoftKey License	8
Add the Sensor	10
GETTING STARTED WITH COLORFLEX EZ	13
Standardization	13
Sensor > Configure Setups	14
Upload	14
Download	14
Edit.....	15
Reading Samples and Standards	15
Product Standard	16
Sample	16
Sensor > Import Logged Reads	16
Retrieve Data.....	16
Sort By Log.....	17
Sort by Setup	17
Select All	17
Copy to Job.....	17
Copy to Database	17
COLORFLEX EZ MAINTENANCE AND TESTING	19
Cleaning of the ColorFlex EZ	19
Diagnostic Test #1: Running Short Term Repeatability	19
Diagnostic Test #2: Running Long Term Repeatability	19
Replacing the Lamp.....	19
COLORFLEX EZ SPECIFICATIONS	21
Operating Conditions	21
Physical Characteristics	21

Conditions of Illumination and Viewing	21
Instrument Performance	22
Regulatory Notice	22
ColorFlex Options and Sample Devices	23
Sample Clamp and Port-Forward Stand (D02-1010-459)	23
Skein Holder Option (02-7396-00)	24
Port Inserts	24
Glass Port Insert	25
UV Port Insert 420nm.....	25
Sample Cup with Ring and Disk Set for Translucent Liquids (CFLX-SC ASSY & LSXE-SC-ASSY)	26
Sample Cup Port Plate.....	27
64mm (2.5in) Glass Sample Cup	27
Sample Cup Opaque Cover.....	27
Ring and Disk Set	27
INSTRUMENT REPLACEMENT, REPAIR, PROBLEMS, AND QUESTIONS	29
Warranty	29
Shipping Claims.....	29
Breakage or Damage	29
Shortage	30
Incorrect Shipment.....	30
Returns.....	31
Packing and Shipping Instruments for Repair	31
When You Need Assistance	32
INDEX	35

ColorFlex EZ Features

The ColorFlex EZ spectrophotometer is a versatile color measurement instrument that can be used for reflectance measurement of products in industries such as paint and textiles. Although the instrument is AC-powered, the small footprint of the ColorFlex EZ allows for portable operation.



Figure 1. ColorFlex EZ

The instrument uses a xenon flash lamp to illuminate the sample. The light reflected from the sample is then separated into its component wavelengths through a dispersion grating. The relative intensities of the light at different wavelengths along the visible spectrum (400-700 nm) are then analyzed to produce numeric results. This is an objective means of quantifying a sample's color.

ColorFlex EZ is available in one geometry – 45° illumination/0° viewing. The label on the back of the instrument provides information on the serial number.

The ColorFlex EZ may be operated using the keypad and display on the instrument or operated while connected to a computer running HunterLab EasyMatch QC software. Having purchased both a ColorFlex EZ and EasyMatch QC software, you have two sources of information on the instrument in addition to this User's Manual: the ColorFlex EZ User's Manual which describes stand-alone operation, and the EasyMatch QC help file, which describes operation of the ColorFlex EZ using the software. Refer to those information sources as required.

ColorFlex EZ Accessories

The following accessories are included with the ColorFlex EZ system and can be found in the provided carrying case:

- **Black glass** - placed at the sample port during standardization of 45°/0° instruments.
- **White calibrated tile** - placed at the sample port during standardization.
- **Diagnostic Green tile** - used to check instrument performance.
- **Standards care card** - gives instructions on how to clean the standards.
- **Certificate of traceability** for the calibrated white tile.
- **USB computer interface cable.**
- **AC adapter, 2 A/9V.**
- **ColorFlex EZ User's Manual**
- **USB Flash disk for Datalog Export**

ColorFlex EZ Installation

The ColorFlex EZ is simple to set up and attach to your computer. Before operating the ColorFlex EZ with EasyMatch QC, you need to install the batteries and connect the instrument to your computer. These steps are outlined below.

1. Unpack the carrying case and remove wrappings and cable ties. Inspect for damage and notify the carrier and HunterLab immediately if any is discovered. Save the packing materials in case it becomes necessary to return the instrument to the factory.
2. Plug the USB cable into the USB port on back left of the ColorFlex EZ.
3. Plug the flat end of the USB cable into the appropriate USB port on the computer.

Install EasyMatch QC Software

Complete the following steps:

1. Log into the system using an account that has **ADMINISTRATOR** privileges for the PC — network or local.
2. Insert the installation CD into the CD-ROM drive. If the system is setup to automatically run CD programs, the menu will appear, and you may skip to Step 5. Otherwise, continue with Step 3.
3. Select the Easy Match QC Icon or from Windows, go to **START > RUN > EZMQC_MENU** and **OPEN**. The following screen will be shown.




Figure 2. EasyMatch QC Installation

4. Select **INSTALL EASYMATCH QC SOFTWARE** and follow the screen prompts.
5. Select **SOFTKEY LICENSE** as the type of key to use with the software.




Figure 3. Software Key License

6. When the EasyMatch QC installation is finished, select the **OPTION BUTTON** next to **YES, I WANT TO RESTART MY COMPUTER NOW** and then **FINISH** to restart the computer and log back in.




Figure 4. Completed Install

7. The CD can now be removed.

Activate the SoftKey License

1. From the Desktop, select the EasyMatch QC Icon or from the Windows Start menu, choose the following to open the software:

START > PROGRAMS > HUNTERLAB > EASYMATCH QC

2. A warning message to activate the license will be displayed as shown below.

Note: EasyMatch QC functions are unavailable before key activation.




Figure 5. No License Warning

3. The SoftKey License is uniquely associated with the sensor serial number and is provided on a thumb drive supplied with EasyMatch QC or via email from HunterLab.
4. Go to **HELP > LICENSE REGISTRATION > ACTIVATION**.
5. Select **ACTIVATE LICENSE**.




Figure 6. Activate License

i. **Option #1: Key ID.**

This method is for copying the ID from an email or writing down the 32-digit code. This requires an internet connection.

- a. From the **CHOOSE AN OPTION** page (Figure 5), select **KEY ID**.
- b. Paste-in or type-in the License Key ID and click **ACTIVATE**.
- c. An acknowledgement will be displayed showing the activation status.

ii. **Option #2: Key File (.skf)**

This method is for using the SoftKey License (.skf file) on the thumb drive.

- a. Place the thumb drive with the SoftKey License in the USB port.
- b. From the **CHOOSE AN OPTION** page (Figure 6), select **KEY FILE (.skf)**.
- c. Browse the USB to find the SoftKey License (.skf file), then click **ACTIVATE**.
- d. An acknowledgement will be displayed showing the activation status.

iii. **Option #3: Sentinel Key**

- a. If the user has a HunterLab USB hardware key, then it can be used with a new sensor on the same computer. Return to Install, Step 5 (Figure 3) and select the **Sentinel Key** to continue.

iv. **Option #4: 30-day trial**

- a. Fill out the registration form provided for the 30-day trial. Connect to the internet. HunterLab will approve the trial and email the SoftKey license back. Follow the directions for Option #1 or #2 to complete.

Figure 7. Request 30-day Trial

Add the Sensor

1. Upon initial startup, the following message will be displayed: **Sensor not yet installed. Please install a sensor to take measurements.** This message will remain until you proceed to the **SENSOR MENU > INSTALL > CONFIGURE** command and install a new sensor.
2. The Sensor Manager appears first:

Figure 8. Sensor Manager

3. Select **ADD SENSOR** to install a new sensor. The Setup Sensor screen allows selection of the instrument model and the communications port. Select **NEXT** when ready.




Figure 9. Setup Sensor

4. Connect the instrument using the USB communications port on the instrument to the computer with EasyMatch QC.
5. Remove the tape covering the reflectance port.
6. Place the desired port plate at the reflectance port and snap it into place.
7. Turn on the ColorFlex EZ by pressing and holding the **READ** Button. Allow the instrument to warm up for two hours prior to standardizing and making measurements.

Getting Started with ColorFlex EZ

Standardization

Standardization on a ColorFlex EZ model with 45°/0° geometry requires reading of the black glass and the calibrated white tile that are contained in the standards box. Standardization can be done through EasyMatch QC (by selecting **SENSOR MENU > STANDARDIZE** or by clicking the **STANDARDIZE** button on the default toolbar) or directly through the ColorFlex EZ firmware.

Standardization sets the top and bottom-of-scale for the neutral axis. During standardization, the bottom-of-scale (zero) is set first. This is done by placing the Black Glass at the sample port. The top-of-scale is then set by reading the calibrated White Tile. Messages on the display prompt through the standardization process as described below.

It is recommended that the ColorFlex EZ be standardized every four to eight hours. The instrument will automatically prompt for standardization after the Standardization Interval set in **SENSOR > SET INTERVAL** has elapsed, so it is a good idea to enter 4 to 8 hours as the Standardization Interval.

Also standardize the ColorFlex EZ any time there is a significant change (greater than 5°F) in ambient temperature. For example, if the ColorFlex EZ is moved an air-conditioned office to an outdoor site that is 90°F, standardization should be repeated once the instrument has had a chance to stabilize under the new temperature. Also standardize the instrument any time there is a change in the optical path, such as changing from using the standard port to the sample cup port.

Also, it is important that the standards used in standardization be treated carefully. Tiles must be kept clean and in good condition. If there is any doubt about the standards being clean, clean them as described in the **Maintaining and Testing ColorFlex EZ** chapter.

Standardization of a ColorFlex EZ is performed as follows:

1. Remove the calibration tiles from the standards box. Check that the tiles are clean and free of dust and scratches. If they are dirty (including marked with fingerprints), clean them as described in **Maintaining and Testing ColorFlex EZ**.
2. **SELECT SENSOR > STANDARDIZE** from the EasyMatch QC. The Standardization screen appears, prompting for the Black Glass. Place the Black Glass on the measuring port with the shiny side of the glass against the port. Check that the Black Glass is flat and in solid contact against the port and that the white dot is facing forward towards the operator.




Figure 10. Place Black Glass at Port

- The ColorFlex EZ reads the **BLACK GLASS** and sets the bottom-of-scale. When it is finished, the screen prompts for the calibrated **WHITE TILE**. Replace the Black Glass with the calibrated White Tile, which is contained in the standards box.




Figure 11. White Tile Reading

- The ColorFlex EZ reads the White Tile and sets the top of scale. When it is finished, the screen indicates that the instrument has been successfully standardized.

Sensor > Configure Setups

For Product Measurements, the Setup Groups window appears to create, edit, or delete setup groups and upload, modify, and download individual setups.




Figure 12. Configure Setups

Upload

When you click **UPLOAD**, the setups stored in the instrument are brought into EasyMatch QC. These setups can then be edited, if desired. The setups can also be saved as a setup group by typing a name into the white box or selecting a name from the drop-down list and clicking **SAVE**. Later, you can select the desired group and click **GET** to retrieve that saved group of setups or **DELETE** to delete the group of setups.

Download

When you click **DOWNLOAD**, the setups in the current EasyMatch QC setup group are sent back to the instrument.

Edit

When you click **EDIT**, the Product Setup Configuration screen appears.

Figure 13. Product Setup

Use the scroll bar next to Setup Number or type a setup number into the box to choose a setup with which to work. If you wish to begin working with the version of the setup that is already resident in your instrument, click **RETRIEVE**. If you wish to work with the version of the setup that is shown on-screen, do not click **RETRIEVE**.

Alter the setup parameters as desired. The parameters and selections available are the same as those configured through your instrument firmware. If you are using a physical standard, you may click the **READ STANDARD** button to read the standard to be saved with the setup using your instrument or **RECALL STANDARD** to recall a standard from your EasyMatch QC database to be saved with the setup.

When all parameters are as desired for this setup, click **UPDATE SENSOR** to send the setup to your instrument. You may retrieve all the setups from the instrument at once using the **RETRIEVE ALL** button or send all the setups back to the instrument at once using the **UPDATE ALL SETUPS** button. Complete the following steps to take individual readings using the ColorFlex EZ:

Note: These instructions apply when Average in the product setup is set to OFF.

Reading Samples and Standards

From the **MEASUREMENTS MENU**, you may read either a standard or a sample. You may also average data. The functions available through the **MEASUREMENTS MENU** are described in the EasyMatch QC Reference Manual.

It is recommended that the instrument be standardized at least once every four hours. Then you may proceed with sample measurement.

Product Standard

A product standard is an object that represents the ideal target color for the product. This object is the one to which others will be compared and deemed acceptable or unacceptable. The product standard may either be a physical item that is measured using the ColorFlex EZ or a set of color values that is entered into EasyMatch QC. .

Sample

A sample is an object that will be measured with the ColorFlex EZ and compared to the product standard. The color of the sample is generally like the color of the product standard.

Sensor > Import Logged Reads

The **SENSOR > IMPORT LOGGED READS** command appears when the connected instrument is a ColorFlex EZ. There must be items saved to the instrument datalog in order to use this command.

The Data log screen appears first.




Figure 14. Data Log Screen

Retrieve Data

Copies of the measurements stored in the instrument's datalog can be retrieved to the Data log screen.




Figure 15. Retrieve Data

Once the measurements are shown, the remainder of the buttons at the bottom of the screen become active and serve the following functions:

Sort By Log

Causes the measurements shown in the Data log window to be sorted by their instrument datalog ID numbers for easier viewing.

Sort by Setup

Causes the measurements shown in the Data log window to be sorted by their instrument setup numbers for easier viewing.

Select All

Selects (highlights) all the measurements shown on the Data log screen so that the **COPY TO JOB** or **COPY TO DATABASE** command may be applied to all of them at once. (The alternative to **SELECT ALL** is to use the Windows convention of clicking on a single item to select it, **CTRL +** clicking on multiple items to select them, or **SHIFT +** clicking on the first and last of a list of items to select all the items between.)

Copy to Job

Places the measurements that are currently selected into the active job. Standards are automatically placed in the job as standards. You are prompted to indicate the standard under which samples should be saved.

Copy to Database

Places the measurements that are currently selected into the database. Standards are saved as standards and samples as samples, and the datalog ID number is used as the item ID.

The measurements copied to the job or EasyMatch QC database are NOT deleted from the instrument's datalog. If you wish to delete them, you must do so manually through the instrument firmware.

ColorFlex EZ Maintenance and Testing

The ColorFlex EZ requires minimal maintenance. This chapter describes cleaning of the instrument and tiles and running diagnostic tests as part of normal instrument maintenance.

Cleaning of the ColorFlex EZ

Clean the outside surfaces of the ColorFlex EZ using a soft cloth. Do not spray liquids directly on the instrument.

Diagnostic Test #1: Running Short Term Repeatability

Short Term Repeatability of your instrument may be tested using the calibrated white tile as follows:

1. Turn the ColorFlex EZ on and allow it to warm up for 2 hours. Meanwhile, clean the calibrated white tile as described and allow the tile to return to room temperature.
2. Follow the instructions given in the ***SENSOR MENU > DIAGNOSTICS*** section to run the short repeatability test that is built into EasyMatch QC.

Diagnostic Test #2: Running Long Term Repeatability

Long Term Repeatability is measured using the Green Tile as follows:

1. If needed, clean the diagnostic green tile and allow the tile to return to room temperature.
2. Follow the instructions given in the ***SENSOR MENU > DIAGNOSTICS*** section to run the long-term repeatability test that is built into EasyMatch QC.

Replacing the Lamp

Lamp replacement requires a trained technician. Contact HunterLab Technical Support to arrange for lamp replacement. Please read “When You Need Assistance” prior to contacting HunterLab.

ColorFlex EZ Specifications

The specifications and characteristics of your instrument are given in this section. For best performance, your instrument should be placed where there is ample work space with medium or subdued illumination and no drafts. The operating conditions (temperature and humidity ranges) are given in the Operating Conditions section below.

Operating Conditions

ColorFlex EZ can be stored in an area with a temperature range of -20°C to 65°C (-5°F to 150°F) for up to 3 weeks and can be operated under temperature conditions of 10°C to 40°C (50°F to 104°F). For specification-level performance, the recommended temperature range is 21-28°C (70-82°F). It may be operated under relative noncondensing humidity conditions of 10% to 90%. Do not leave ColorFlex EZ in an area where temperature or humidity extremes are possible.

Physical Characteristics

Weight	4.5 kg (9.9 lb.)
Dimensions	Height: 16 cm (6.3 in) Width: 13 cm (5.1 in) Depth: 36 cm (14.2 in)
Communications Interface	3 USB 2.0 ports
RFI Compliance	FCC Class A (Commercial), IEC, or equivalent
Safety Compliance	UL, CSA, IEC, or equivalent
System Power	100 to 240 VAC, 47 to 63 Hz

Conditions of Illumination and Viewing

Light Source	Pulsed xenon
Source UV content	Match to D65 with CIE rating of CC or better
Lamp Life	>1 million flashes
Geometry	Directional annular 45° Illumination/0° Viewing
Detector	Sealed optics; 256-element scanned array & high-resolution concave holographic grating
Port Diameters/Sample View Diameters	31.8 mm (1.25 in)/25 mm (1.0 in)

Instrument Performance

Spectral Data	Range: 400-700 nm Reporting Interval: 10 nm
Bandwidth at Half-height	10 nm
Photometric Range	0-150% reflectance
Measurement Speed (at 23°C)	≤1.5 seconds
Measurement Storage Capacity	2000 spectral readings as sample 250 spectral or tristimulus standards with Pass/Fail Tolerances 100 product setups

Note: Every attempt at accuracy is made, but specifications are subject to change without notice.

Note: Use of this equipment in a manner not specified by the manufacturer may impair the protection afforded by the equipment. Danger of electric shock if liquids are spilled and fire if volatile or flammable liquids are spilled. Use care when measuring liquid samples.

Regulatory Notice

A copy of the Declaration of Conformity according to ISO/IEC Guide 22 and EN 45014 follows:



Figure 16. Sample Clamp and Stand

Skein Holder Option (02-7396-00)

This is a device for measuring yarn skeins. Wind the yarn around the skein holder in multiple taut layers until it is effectively opaque and is as flat as possible.

Place the skein holder at the sample port and back it with the sample clamp or a white backing tile to provide a consistent background and pressure. Make several measurements of the skein, rotating the holder 90° between measurements and averaging the readings for the result.



Figure 17. Skein Holder

Port Inserts

Special port plate inserts with various-sized openings are available for use with samples of differing sizes. Some inserts are also available with glass covers to protect the inside of the instrument from sample and airborne particles.



Figure 18. Port Inserts

Installation:

Remove the standard port insert by twisting it until it pops up. Replace it with the desired port insert. Make sure that the flat portion of the insert is facing outward with the magnetized portion toward the interior of the sensor.

Sample Measurements:

Standardize the instrument. Any change in port size must be followed by standardization on both the Black Glass and the Instrument White tile. When using a glass port insert, standardize the instrument with the glass insert in place. However, you should check measurements of the green and white tile and perform other diagnostics *without* the glass insert. Difference measurements using a glass port insert will be more accurate than absolute measurements.

Glass Port Insert

This option provides a snap-in port plate with an opening one inch in diameter that is covered with glass. The glass insert protects the inside of the instrument from sample and airborne particles.

Installation:

Remove the standard port insert by pulling it out of the port insert retainer. Replace it with the glass port insert. Make sure that the flat portion of the insert is facing outward with the beveled portion toward the interior of the sensor.

Sample Measurements:

Standardize the instrument with the glass in place. Check the measurements of the green and white tiles before performing any other diagnostics without the glass cover. Difference measurements using a glass cover will be more accurate than absolute measurements.



Figure 19. Glass Sample Port

UV Port Insert 420nm

This port insert contains a 420-nm filter to prevent light below 420 nm from hitting the sample. This is important when samples have been treated with UV enhancing chemicals.

Installation:

Remove the standard port insert by pulling it out of the port insert retainer. Replace it with the UV port insert. Make sure that the flat portion of the insert is facing outward with the beveled portion toward the interior of the sensor.

Sample Measurement:

Then standardize the instrument with the UV filter in place. Check the measurements of the green and white tiles before performing any other diagnostics without the UV filter.

Sample Cup with Ring and Disk Set for Translucent Liquids (CFLX-SC ASSY & LSXE-SC-ASSY)

This set provides a device for holding powders, pellets, granules, and translucent liquids at the sample port. It includes a 64mm (2.5-inch) glass sample cup, plastic ring and ceramic disk set, special port plate and sample cup cover. The ring and disk set are provided for translucent liquids, transparent liquids, and semi-solids where the path length is related to color concentration and must be fixed.



Figure 20. Ring & Disk Set

Installation:

First, the sample cup port plate should be installed to hold the sample cup in the proper position. Place the ring inside the sample cup and fill the cup with sample to a level above the ring. Then place the ceramic disk on top of the sample (until it rests on top of the ring) with the white portion facing the sample. Cover with the sample cup cover to exclude external light.

Sample Measurement:

Measure the sample at the sample port through the glass bottom of the sample cup with the instrument in the port-up orientation.

Parts:

Accessory Parts	ColorFlex EZ	LabScan
Sample Cup Port Plate	04-6622-00	02-1010-316
64mm Glass Sample Cup	04-7209-00	04-7209-00
Sample Cup Opaque Cover	04-4000-00	04-4000-00
Ring and Disk Set	02-7579-00	02-7579-00

Sample Cup Port Plate

This port insert is specially designed to accommodate the 64mm (2.5-inch) glass sample cup and to hold it in the proper position at the measurement port. Remove the standard port insert by pulling it out of the port insert retainer. Replace it with the sample cup port insert, then standardize the instrument.

64mm (2.5in) Glass Sample Cup

The sample cup is ideal for the analysis of powders, granules, pellets, and translucent samples. It is recommended for use with the instrument in the port-up orientation. The cup can be filled with the desired sample and placed at the measurement port so that the sample is measured through the glass bottom of the cup.

Sample Cup Opaque Cover

The opaque cover provides a light trap to exclude the interference of external light on a sample contained in the 2.5-inch sample cup. It is placed over the filled cup at the sample port as shown below.



Figure 21. Opaque Cover

Ring and Disk Set

The ring and disk set includes a plastic ring and ceramic disk designed for use inside the 2.5-inch glass sample cup (sold separately). The ring and disk set are recommended for use when measuring translucent liquids, transparent liquids, and semi-solids to exclude external light and provide a consistent white background and sample path length.

Place the ring inside the sample cup and fill the cup with the sample to a level above the ring. Then place the ceramic disk on top of the sample (until it rests on top of the ring) with the white portion facing the sample. If desired, cover the entire assembly with the opaque cover (sold separately). Measure the sample at the sample port through the glass bottom of the sample cup with the instrument in the port-up orientation.



Figure 22. Ring and Disk Set

Instrument Replacement, Repair, Problems, and Questions

The following HunterLab policies are described in this chapter:

- Warranty
- Claims
- Returns/Service
- Technical Assistance.

Warranty

HunterLab warrants that all instruments it manufactures are free from defects in material and workmanship under normal use. This warranty is limited to repairing or replacing any defective hardware or software that may cause the instrument to perform outside of its specified tolerances. This warranty is one year from date of shipment of new instruments and two months from the date of shipment of repaired instruments.

Note that printers and computers are covered under the original manufacturer's warranty.

The warranty is void if the user has made unauthorized repairs, improperly installed, operated, or subjected the instrument to conditions outside of the specifications in the product documentation.

The HunterLab warranty does not cover consumable items such as lamps, fuses, batteries, etc.

Questions concerning operation, maintenance, or repair of your equipment can be directed to the Service Department at Service@hunterlab.com. Additional information can be obtained at <http://support.hunterlab.com>.

Shipping Claims

All materials are sold F.O.B. from Reston, Virginia (unless otherwise specified) and HunterLab responsibility ends upon delivery to the first carrier. All claims for loss or damage must be rendered by the consignee against the carrier within fifteen days of receipt of goods. A copy of this notice must also be forwarded to HunterLab within five days of its receipt.

Breakage or Damage

According to the contract terms and conditions of the carrier, the responsibility of the shipper ends at the time and place of shipment. The carrier then assumes full responsibility. Perform the following procedures if your instrument arrives broken or damaged.

Freight or Express

1. Notify your local carrier.
2. Hold the damaged goods with their container and packaging for inspection by the examining agent. Do not return any goods to HunterLab prior to inspection and authorization of the carrier.
3. File a claim against the carrier. Substantiate this claim with the examining agent's report. A certified copy of our invoice is available upon request. The original B/L is attached to our original invoice. If the shipment is prepaid, write for a receipted transportation bill.
4. Advise HunterLab regarding replacement.

Parcel Post Shipment

1. Notify HunterLab at once in writing, giving details of the loss or damage. This information is required for filing a claim.
2. Hold the damaged goods with their container and packaging for possible inspection by postal authorities.
3. Advise HunterLab regarding replacement.

United Parcel Service

1. Contact your local UPS office regarding damage and insurance claim. Each UPS office has a different method of handling these occurrences and yours will advise you of its procedures.
2. Retain the container and packaging.
3. Notify HunterLab at once for replacement.

Shortage

Perform the following procedure if your order appears to be missing items.

1. Check the packing list notations. The apparent shortage may be a back ordered item and may be marked as an intentional short-ship.
2. Re-inspect the container and packing material, particularly to locate smaller items.
3. Ascertain that the item was not removed by unauthorized personnel prior to complete unpacking and checking.
4. Notify HunterLab immediately of the shortage in writing.

Incorrect Shipment

Perform the following procedure if material received does not correspond with your order.

1. Notify HunterLab immediately, referencing order number and item.
2. Hold incorrect items until return shipping instructions are received.

Returns

A service request order (SRO) number is required before any items can be returned to HunterLab. Contact HunterLab's [Order Processing Department](#) to obtain an SRO for damaged or incorrect parts, or [the HunterLab Service Department](#) to obtain an SRO to return an instrument for service.

Do not return any damaged or incorrect items to HunterLab until all shipping instructions are received.

Note: HunterLab must be notified within fifteen days or we cannot accept responsibility for damaged or incorrect items.

HunterLab offers complete repair service for all instruments it manufactures. Call HunterLab for the service facility nearest your location. If your equipment is not functioning properly, contact the HunterLab Service Department for maintenance or repair instructions. Many times, this on-the-spot diagnosis is all that is required.

If repair is required, HunterLab offers two means of servicing. Instruments may be returned to a HunterLab service facility for repair or a HunterLab Service Department technician can come to your location to perform on-site repair. For schedule and terms for on-site repairs by trained service technicians, call the HunterLab Service Department. Please read "When You Need Assistance" prior to contacting HunterLab.

The customer is responsible for incoming and outgoing freight charges for instruments being returned to HunterLab for all repairs, including warranty repairs.

Packing and Shipping Instruments for Repair

Please regard the following instructions when packing your instrument to return it to HunterLab for repair. **Proper packing is crucial.** These instructions do not replace the recommended professional packaging for your instrument, but may assist in eliminating the need for a shipment claim due to faulty packaging. Purchasing freight insurance does not guarantee a successful damaged shipment claim if the carrier determines the instrument was not packaged properly.

- All instrument tiles, the didymium filter (if included), black glass or light trap, power supply, power cords, and cables for the instrument should be included in your shipment. **Your repair estimate will be delayed if the instrument tiles are shipped separately later.**
- Remove the sample clamp (if you have one) from the instrument before packing.
- Cover the measurement port. If applicable, also cover the transmission port and tape the transmission compartment door closed. **Do not use duct tape.** Painter's tape is preferred, as it will not leave residue on the instrument.
- Insert the instrument into an anti-static or plastic bag prior to placing it in the carton. The bag will aid in keeping packing material out of the instrument.
- Place the bag-wrapped instrument into a new carton which includes, at a minimum, **6 inches of packing material** (preferably foam) around the instrument. Styrofoam peanuts should not be used as packing material for instruments, as they can suspend items weighing only up to 5 pounds. Observe the information listed on the bottom of most cartons with regard to burst strength and gross weight limits. Single wall cardboard cartons should not

be used. (A proper packing carton with packing material may be purchased from HunterLab, if desired.)

- Insure the shipment.
- Provide an itemized packing list of all contents of the shipment.
- Label the carton(s) as follows:

Hunter Associates Laboratory Inc.
Attn: SRO # _____
11491 Sunset Hills Road
Reston, VA 20190
U.S.A.

When You Need Assistance

When you have a problem with an instrument or software or need technical advice concerning a specific application, you may consult the support website (support.hunterlab.com). There are numerous articles on applications, operations, instrument accessories, troubleshooting and more. This is available 24/7. If you don't find the information that you require, then you can open a support request on the website. Please include the following information when corresponding with HunterLab.

1. The type of sensor in use.
2. The serial number of the instrument (usually found on a tag on the back or bottom of the sensor, or inside the transmission compartment).
3. The type of software you use to access the sensor output (EasyMatch QC), the version of the software (seen after choosing **HELP > ABOUT**), the operating system, and the brand and type of computer.
4. The specific nature of the problem, including the exact error message received or the number of units the sensor reads **OFF** from the standard tiles.
5. The steps performed prior to the start of the problem.
6. Steps already performed to reconcile the problem and/or results of any diagnostics.
7. The type of product being measured.
8. Operating environmental conditions under which the instrument is normally used, such as temperature, humidity, dust, fumes, etc.
9. Whether the instrument has recently been moved or the computer reconfigured.
10. The name(s) of any HunterLab personnel with whom you have previously discussed the problem.

To place an order, for prices on instruments, software, or replacement parts, or to return damaged or incorrect parts, ask for the Order Processing Department. For applications advice or for help in correcting instrument or software problems, ask for Technical Support. To return instruments to HunterLab for service, or to ask questions about the servicing or recalibration of instruments, ask for the HunterLab Service Department. To speak with HunterLab, please call 703-471-6870.

The mailing address for HunterLab headquarters is given below. Customers outside the United States should contact their HunterLab distributor for initial assistance.

Hunter Associates Laboratory, Inc.
11491 Sunset Hills Road
Reston, Virginia 20190 U.S.A.

Index

- AC adapter, 6
- Accessories, 6
- Black glass, 6
- Broken instruments, 29
- Certificate of traceability, 6
- Claims, 29
- Cleaning the ColorFlex EZ, 19
- Damaged instruments, 29
- Features, 5
- Green tile, 6
- Illumination, 21
- Incorrect shipment, 30
- Installation, 7
- Instrument Performance, 22
- Lamp Replacement, 19
- Maintenance, 19
- Operating Conditions, 21
- Options, 23
- Order processing department, 32
- Packing instrument, 31
- Physical characteristics, 21
- Problems, 28
- Questions, 28
- Reading, 15
- Regulatory notice, 22
- Repair, 28
- Repeatability test, 19
- Replacement, 28
- Returns, 31
- Sample, 16
- Sample Clamp & Stand, 23
- Sample Cup with Ring and Disk Set, 26
- Sample devices, 23
- Service Department, 33
- Shipping claims, 29
- Shipping instrument, 31
- Shortage, 30
- Skein Holder, 24
- Specifications, 21
- Standard, 16
- Standards care card, 6
- Technical Support, 33
- Testing, 19
- USB computer interface cable, 6
- Viewing, 21
- Warranty, 28
- When you need assistance, 32
- White calibrated tile, 6