

User's Manual for Aeros[®] with EasyMatch[®] Essentials



Hunter Associates Laboratory
11491 Sunset Hills Road
Reston, Virginia 20190 USA
www.hunterlab.com

A60-1018-193 Version 2.0
For EasyMatch[®] Essentials Version 1.05.0090 and Above

Copyrights and Trademarks

This documentation contains proprietary information of Hunter Associates Laboratory, Inc. Its reproduction, in whole or in part, without express written consent of Hunter Associates Laboratory, Inc. is prohibited.

Aeros and EasyMatch are registered trademarks of Hunter Associates Laboratory, Inc.

Windows is a registered trademark of Microsoft Corporation in the United States and other countries.



Caution: If the equipment is used in a manner not specified by the HunterLab, the overall safety may be impaired. - The instrument is for indoor use only and not suitable for a wet location.



Caution: There is a potential of a UV Light hazard in using this instrument. Please avoid looking directly at the light. The frequency of this flashing light is in the range of sensitivity for those prone to epileptic seizures. When reading a sample, the illuminated spot flashes in the range of 5 Hz to 6 Hz. User discretion is advised.



Safety Notes

For your safety when using the Aeros, you should pay attention to the following types of statements in this User's Manual:

- General safety instruction that should be observed at all times while operating the instrument.
- Specific safety instruction critical to the type of instrument operation being explained in the manual where the caution appears.
- Use of this equipment in a manner not specified by the manufacturer may impair the protection afforded by the equipment.
- Danger of electric shock if liquids are spilled and fire if volatile or flammable liquids are spilled. Use care when measuring liquid samples.
- This instrument sensor moves up and down during standardization and creating a measurement profile. Please keep fingers and other items out of the way of the sensor.
- The turntable will rotate if turned on in Read Options > Measurement configuration. Please take care to remove fingers, jewelry and clothing to prevent damage.
- The Aeros is for indoor use only at an altitude of up to 2000m and pollution degree 2.

Contents

SETTING UP THE AEROS	7
Standard Accessories	7
Power Jack.....	7
Power Switch.....	8
Keyboard And Mouse.....	8
Front And Rear USB Connectors	8
Ethernet Port.....	8
Moving The Unit.....	9
TAKING A SIMPLE MEASUREMENT.....	11
What Is HunterLab Aeros & EasyMatch Essentials?	11
Connecting The Sensor And Taking A Measurement.....	11
Standardize	11
Create a New Workspace	13
Read Sample	15
NAVIGATING THE ESSENTIALS SCREEN.....	17
Tools: Read.....	17
Tools: View Flippers	17
Tools: Status Bar.....	18
INFORMATION	18
STANDARDIZATION STATUS	18
JOBS STATUS.....	18
WORKSPACE STATUS	18
Tools: Recall Measurements.....	18
Tools: Move Sensor To Top.....	19
Tools: View Options	19
Tools: Workspace & System Settings.....	19
Tools: Jobs.....	19
TOOL BAR: JOBS FUNCTION	21
Jobs	21
JOBS > NEW	21
JOBS > OPEN	21
JOBS > SAVE & SAVE AS	22
JOBS > PRINT.....	22
JOBS > DATA MANAGEMENT.....	23
JOBS > HELP	29
JOBS > ABOUT.....	29
TOOL BAR: WORKSPACE & SYSTEM SETTINGS.....	33
Workspace > Color Scales	33
Workspace > Read Options.....	37
READ OPTIONS > OPTIONS	37
READ OPTIONS > MEASUREMENT CONFIGURATION	40
Workspace > Standard And Tolerances	42

AUTOTOLERANCES SETTING	45
TOLERANCES	46
Workspace > Views	46
Workspace > New Workspace	47
System Settings: Standardization.....	47
System Settings > Diagnostics.....	49
TESTING THE AEROS FOR COLORIMETRIC REPEATABILITY.....	49
READING THE GREEN TILE	52
RUN AUTO DIAGNOSTICS	52
VALIDATE	53
ADVANCED TESTS	53
PREDICTIVE TESTS.....	57
System Settings > Preferences.....	59
PREFERENCES > GENERAL.....	59
PREFERENCES > PRINT	61
System Settings > User Manager	62
TOOLBAR: OPTIONS (VIEWS).....	65
Views > EZ View	65
Views > Color Data Table	66
Views > Spectral Data Table.....	68
Views > Spectral Plot.....	69
Views > Trend Plot	70
Views > Color Plot	71
RECALL STANDARD	73
SPECIAL FUNCTIONS	75
Auto-Exporting Data Through A Network Connection	75
CONNECT TO A NETWORK HUB USING AN ETHERNET CABLE	75
Auto-Exporting Data Via Direct Connection Between Aeros And A Computer	77
Tips & Tricks: Recover Unsaved Measurement Data.....	81
HunterLab File Service Package	81
To Install the File Service Package.....	81
Setting Up File Storage From Aeros (Client) Side	82
SPECIFICATIONS	85
Operating Conditions	85
Physical Characteristics	85
Conditions Of Illumination And Viewing.....	86
Instrument Performance.....	86
Measurement.....	87
Regulatory Notice	88
FEATURES, ACCESSORIES & MAINTENANCE	89
Aeros Maintenance & Safety	89
CLEANING THE AEROS	89
CLEANING THE WHITE TILE.....	89
CLEANING THE BLACK GLASS AND GREEN TILE	89

User's Manual for Aeros and EasyMatch Essentials v 2.0

POWER REQUIRED	89
SAFETY	90
WHEN YOU NEED ASSISTANCE	91
INDEX	93
TABLE OF FIGURES	95

Setting Up the Aeros

The Aeros system contains two major components – the sensor head and the turntable base. The sensor head contains the spectrophotometer with an LED light source, distance measuring components and a touch screen for Essentials. The turntable base provides the electronics and the mount for the sensor with automated vertical positioning and horizontal manual adjustment. The turntable rotates the product under the optical sensor.

Standard Accessories

- Calibration Box with calibrated white tile, black glass and green diagnostics tile
- Certificate of Traceability
- Power Supply
- Aeros Quick Start Guide
- 12-in and 6-in Sample Dishes
- USB Flash Drive



Figure 1. Calibration Box

Power Jack

The instrument is supplied with a 24 VDC (3.75A) power supply. The power supply is plugged into the back of the instrument as shown along with the Ethernet port and the USB port.




Figure 2. Rear View with Connectors

CAUTION: Use only the power cord included with this instrument or a replacement obtained from HunterLab. Be certain that the power cord is in good condition before connecting it.

Power Switch

To turn the instrument on, press the rocker switch on the back of the instrument.

Keyboard And Mouse

The Aeros works with the following keyboard and mouse:

- L02-1017-434 Wireless keyboard and mouse kit.

To use this accessory, turn the power off. Plug in the micro USB adaptor to the rear of the instrument and then attach the nano-receiver for the keyboard into the USB port. Install the batteries into the keyboard/mouse and turn the power back on.

Front And Rear USB Connectors

There are two USB connectors on the Aeros. The one in the front is typically used for exporting jobs and workspaces, backing up the instrument and updating software. The USB port on the back of the instrument is typically used to connect a printer or a keyboard to the Aeros.




Figure 3. USB Port on Rear of Instrument

Ethernet Port

This port is used to connect the Aeros to:

- Computer or to a network with the purpose of sending data (ASCII) to a server.
- Connect with EasyMatch QC and EasyMatch ER
- Remote Support
- Network printer.

Moving The Unit

Use care in moving the instrument which weighs 50 lbs. Training on lifting heaving objects is recommended. To lift from the box, a two-person lift is suggested. When moving it from a table to another location, balance this instrument by lifting from behind and by the base.

If moving the Aeros any distance, please secure the optics using **DIAGNOSTICS > ADVANCED > PARK** for Shipping.

Taking a Simple Measurement

What Is HunterLab Aeros & EasyMatch Essentials?

Aeros is a reflectance-only color measuring instrument capable of measuring the color of irregularly shaped/textured products. All samples are measured by placement in a tray or container under the sensor head. With features like Auto Height Positioning, a rotating sample platform, large touch-screen display, and smart communications, the innovative Aeros can measure products like coffee beans, snack foods, plastic pellets, even industrial minerals.

Connecting The Sensor And Taking A Measurement

After unpacking and setting up the instrument, turn on the power using the rocker switch on the back of the instrument base.


Once inside the software, the main measurement screen is displayed – **COLOR DATA TABLE (D65/10)**.

Figure 4. Measurement Screen



The instrument is automatically connected and this is reported on the status bar. Next, the unit must be Standardized.

STANDARDIZE

- Press the **WORKSPACE** icon  and select **STANDARDIZATION** User will be prompted to remove all samples from under the sensor. .
Alternately, press the **STANDARDIZATION** status area on the status bar to initiate Standardization.
- **MOVE SENSOR TO TOP:** Press **OK** to move the sensor to the highest position.

User's Manual for Aeros and EasyMatch Essentials v 2.0



Figure 5. Sensor Moves to Top

- **READ BLACK GLASS:** When the sensor stays at the top, attach the standardization box to the sensor. Then attach the black glass to standardization box and press **READ**.



Figure 6. Read the Black Glass for Bottom-of-Scale

- **READ WHITE TILE:** Remove the black glass and attach the calibrated white tile to standardization box. Press **READ** to continue.




Figure 7. Read the White Tile for Top-of-Scale.

- Remove the calibration box when standardization is completed. Click **OK**



Figure 8. Completed Standardization

- If the user has selected the **GREEN TILE READING DURING STANDARDIZATION** under **WORKSPACE > PREFERENCES**, the user will be prompted to read the green tile too. The target values are the ones on the back of the green tile.
- Standardization is updated and reported as **STANDARDIZED** on the bottom of the screen.



Figure 9. Standardization Status Bar

CREATE A NEW WORKSPACE

- Create a workspace for your product. To **CREATE A NEW WORKSPACE**: From the **WORKSPACE** Menu, press **NEW WORKSPACE** and enter a name for this Workspace. A new job will be opened together with this new workspace.




Figure 10. New Workspace

Note: The active workspace name is displayed in the lower right corner of the screen.

- **CONFIGURE THE NEW WORKSPACE:** With this new workspace open, you can start to change the settings in Color Scales, Read Options, Standard and Tolerances, Views and View Options. All these changes will be automatically saved in this new workspace.
- **DEFAULT SETTINGS:** Following are the default settings for a new workspace:

Table 1. Workspace Parameters

Parameter	Selection
Color Scales	CIE L*a*b*
Illuminant	D65/10
Indices	None
Differences	None
Read Options > Options	Prompt for Sample Name, Auto Save
Read Options > Measurement Configuration	Turntable ON, Auto Height selected, Measurement time = 5 sec
Standard and Tolerances	None
Views	Color Data Table only
View Options for Color Data Table	Latest Data First selected, Precision = 2


- When **AUTO HEIGHT** is selected as default, the sensor head will adjust its position automatically for each sample measurement.
- Users can also choose to read samples with the sensor head at a fixed distance. **FIXED HEIGHT:** No matter what height the sensor is detected, it will use the height that customer enter here to do color calculation. This is used when the sample is hard to detect the height (e.g. sample has absorption at 880nm. The laser of height sensor is absorbed by sample)
- **CREATE SAMPLE PROFILE** is covered in **WORKSPACE > READ OPTIONS > MEASUREMENT CONFIGURATION** portion of this manual. (Create sample profile is to use the detected height for all sample measurements. If sample height out of the desired range (+/- 0.5inches), then it will prompt the error. Sensor will use the detected height to calculate each sample measurement.



Figure 11. Read Options: Measurement Configuration

- Now your instrument is ready to read your product under this new workspace. If you'd like to start a new job for this product, you can press **NEW JOB** and load this configured workspace to continue.

READ SAMPLE

- **PREPARE SAMPLE:** Place the sample dish on the turntable.
- **READ SAMPLE:** To read a sample, press the green lightning bolt. 
- **SAMPLE NAME:** The default sample name is Sample + numerical increment. To customize the sample name, go to **WORKSPACE > READ OPTIONS > PROMPT FOR SAMPLE/STANDARD NAME**. Select the **PROMPT for SAMPLE NAME** to manually input the name during the measurement cycle. Or change the default Sample Name to another name for numerical sequence. Press **APPLY** when done.




Figure 12. Configure Sample Name

Note: The standard/sample name in Essentials should not be empty and should not contain any of the characters , ; ' " + = ? * < > \ /.

- **MAIN MEASUREMENT SCREEN:** The Color Data Table view shows the configured Color Scale results for the standard and sample measurements in the job. The configured tolerances can be applied to the Job and Pass/Fail results will also be displayed. To add the color differences, indices and tolerances to the Color Data Screen, see **WORKSPACE > COLOR SCALES** and **WORKSPACE > STANDARD AND TOLERANCES**.
- To add a product standard and tolerances, see **STANDARDS AND TOLERANCES**. To change the color scale, etc., see **WORKSPACE > COLOR SCALES**. The setups are saved automatically in Workspace.
- A long press on the Sample name will show a menu with the following options:
 - **SET AS STANDARD** – to set the sample as Standard.
 - **RENAME** – to rename the sample.

- **DELETE** – to delete the sample.



Figure 13. Changing, Renaming or Deleting a Sample

- A long press on the Standard name will show a menu with the following options:
 - **EDIT** – to edit the standard. If Edit is selected, the Workspace > Standard and Tolerances dialog box is presented to allow for editing the name, assigning tolerances or changing the type of standard.
 - **DELETE** – to delete the standard. The deleted Standard is reverted into the samples list with its original name.



Figure 14. Edit/Delete a Standard

Navigating the Essentials Screen

The EasyMatch Essentials Tools and Status features are shown below.




Figure 15. User Interface Screen for Aeros & Essentials

Tools: Read



Samples are read using this key. This tool can be moved around the screen by pressing and moving the icon.

Tools: View Flippers



Switching between Views can be accomplished by using the semi-transparent **NEXT** and **PREV** buttons placed at the side edges of the screen or by swiping left or right with two fingers on the screen.

Tools: Status Bar

INFORMATION



The sensor type and serial number is shown at the bottom left side of the System Bar when the *i* is pressed. When application security is enabled and the user logs into Essentials, the User Account will also be shown in the Information box.



Figure 16. Sensor Serial Number

STANDARDIZATION STATUS

The current state of standardization is reported. To initiate standardization, one can press on the **SENSOR STATUS** to open the standardization dialog.

JOBS STATUS

JOB Status is reported on the bottom right side of the System Bar. To open a Job, one can click on the **JOB NAME** in the status bar.

WORKSPACE STATUS

WORKSPACE Status is reported on the bottom right side of the System Bar. . To load a Workspace, click on the **WORKSPACE NAME** in the status bar.

Tools: Recall Measurements



Allows for quick recall of Standards and Samples.

Tools: Move Sensor To Top



This tool moves the sensor head to the top.

Tools: View Options



This menu shows the configuration options for the active view. A total of six views are available. Each view shows a different option. Views can be added or removed in **WORKSPACE > VIEWS**.

Tools: Workspace & System Settings



The **WORKSPACE** menu sets up the data screen with **MEASUREMENT COLOR SCALES, READ OPTIONS, STANDARDS, TOLERANCES** and **VIEWS**.

SYSTEMS SETTINGS initiates **STANDARDIZATION, DIAGNOSTICS, PREFERENCES**, and the **USER MANAGER** for System Security.

Tools: Jobs



A **JOB** is a collection of all the sample measurements and a workspace used for a task, product, or customer. Jobs are the readings of EasyMatch Essentials. Jobs can be created for many different reasons, such as to hold data for a certain customer or a specific product line. Each operator may maintain their own job with preferences or create separate jobs for different operations.

A **WORKSPACE** is a collection of the measurement parameters for a job along with tolerances and the standard, i.e. analogous to word processing documents containing text and formatting. Each job has only one workspace.

Tool Bar: Jobs Function

Jobs Icon



Under the Job function, the following tasks can be accomplished:




Figure 17. Job Menu

Jobs

JOBS VS. WORKSPACE: A job consists of samples measured according to a specific workspace.

A workspace is a template with measurement conditions such as **STANDARD & TOLERANCES, COLOR SCALE, INDEX, ILLUMINANT**, etc. There can be only one workspace associated with a job. The main tool bar provides the options to create a **NEW JOB, OPEN AN EXISTING JOB** and **SAVE A JOB**.

JOBS > NEW

When click **NEW JOB**, the **LOAD WORKSPACE** dialog will pop out. The default selected workspace is the current workspace. User can change the workspace and click **LOAD** then the selected workspace is opened in the new job. Once the Workspace has been loaded to the New Job, the Workspace name associated with this job cannot be edited.

- **Shortcut:** press Workspace name in the lower right status bar

JOBS > OPEN

OPEN a saved Job: A list of available jobs under the current path are displayed for selection. If the job that is needed exists in another folder, then it is an option to change the folder. When the job to be opened is displayed select the file and press **OPEN**. We also provide a short cut for Jobs: Open in main tool bar.

- **Shortcut:** press job name in the status bar on the lower right.

User's Manual for Aeros and EasyMatch Essentials v 2.0




Figure 18. Open A Job

JOBS > SAVE & SAVE AS

SAVE THE JOB under the desired name: To save a job, **SELECT THE FOLDER, NAME THE JOB** and **SAVE THE JOB** contents into a file. These files have an '.ezm' or a CSV extension. There will be a default name filled in Filename box as date&time&instrument#&workspace. You can edit it if needed.




Figure 19. Save A Job

JOBS > PRINT

PRINT AN OPEN JOB using the parameters set up under **WORKSPACE & SYSTEMS SETTINGS > PREFERENCES**.




Figure 20. Workspace & System Settings > Preferences > Print

Drivers included in the Aeros are shown below. Additional printer drivers can be added under **WORKSPACE > DIAGNOSTICS > ADVANCED**.

Table 2: Printer Drivers Included

Printer	Driver
Canon	Canon Print Service 4.4+
HP	HP Print Service Plugin 4.1+
Epson	Epson Print Enabler 4.4+
Konica Minolta	Konica Minolta Print Service Plugin 4.4+
Kyocera	Kyocera Print Service Plugin 4.4+
Lexmark	Lexmark Print Service Plugin 4.4+
Sharp	Sharp Print Service Plugin 4.4+
Xerox	Xerox Print Service Plugin 4.4+

- Printing can be downloaded to a pdf file by selecting, **SAVE AS PDF**. Once this is selected, the parameters for the output are presented. Please save the file to the download folder.




Figure 21. Save as PDF.

Figure 22. Save PDF to Downloads File




Figure 23. Download File Contents

JOBS > DATA MANAGEMENT

Standard(s) and sample measurements are saved in Job files and database along with the sensor information. The saved measurements are also associated with a respective Workspace and Job.

The **DATA MANAGEMENT** contains the features to **RECALL, IMPORT, EXPORT, EMAIL A JOB** and **BACKUP/RESTORE**.

- **RECALL** the measurements from the database.

- **IMPORT** a selected Job(s), Standard(s), Workspace(s)Diagnostics, photos for logo print setup and others from a USB flash drive.
- **EXPORT** the Job(s), Standard(s), Workspace(s), Diagnostics, pdf reports and others to a USB flash drive.
- **EMAIL** the selected Job(s), pdf reports and other files.
- **DELETE** Job(s), Standard(s), Workspace(s), Diagnostics, pdf reports and others.
- **BACKUP** the Hunter Lab folder (all jobs, database and user manager settings) into a USB Flash drive.
- **RESTORE** the Hunter Lab folder (all jobs, database and user manager settings) from a USB Flash drive.




Figure 24. Data Management Menu

Jobs > Data Management > Recall

Recall measurements that have been stored to a job. This feature enables the user to recall the Standard/Sample(s) stored in the Database into the current running Job. A shortcut to this function is available on the toolbar.

Click **RECALL** option in Data Management, a dialog will be displayed where the user can recall the measurements from the database by **SELECTING THE TYPE**:

- Show The Standard Associated With A Specified **CATEGORY**.
When this options is selected, the standards list is filled with the standards associated with the selected Category from the **SELECT STANDARD CATEGORY** list. All matching samples are filled into the sample list. To narrow the list, the user can select the samples only associated with the selected standard.
- Show All The Measurements In The Selected **JOB**
When this option is selected, the standard and samples for that job are listed.
- Show The **STANDARDS/SAMPLES** In The **CURRENT WORKSPACE**
When this option is selected, the list is filled with the standards and samples associated with the current workspace. To narrow the sample list, the user can select the samples only associated with the selected standard.
- Show The **IMPORTED STANDARDS**

When this option is selected, the standards list is only filled with the standards imported into the database. After selecting the measurements click **RECALL** placing them into the active Job.




Figure 25. Recall Measurements

Jobs > Data Management > Import

This feature allows the user to import the below data from a USB flash drive into the instrument. Data can be one file or multiple files. All selected files should be in the same file path location. The following data can be imported: **JOBS, STANDARDS, WORKSPACE, DIAGNOSTICS,** and **OTHERS.**

- **IMPORT JOB**

This option allows the user to browse and select a Job file(s) (.ezm) from the USB flash drive and import into the instrument. If a file name already exists, then the name will be incremented numerically. To use these functions, a USB flash drive must be present in the port.




Figure 26. Import Job

- **IMPORT STANDARD**

This option allows the user to browse and select a Standard(s)(extension .std) from the USB flash drive and import into the database. If required, the Standard Name can be changed.

- **IMPORT WORKSPACE**

This option allows the user to browse and select a Workspace(s) (extension .wsp) from the USB flash drive and import into the database. If the workspace already exists, then the user is prompted to specify a different name.




Figure 27. Import Workspace

- **IMPORT DIAGNOSTICS**

This option allows the user to browse and select a Diagnostics file from the USB flash drive for import into the instrument database.

- **IMPORT OTHERS**

This function is available to import other items such as a logo for the printed report. Once the logo is imported, go to **WORKSPACE & SYSTEM SETTINGS > PREFERENCES > PRINT** to add the logo to a report.

Jobs > Data Management > Export

This feature allows the user to export the below data from instrument into a USB flash drive. Data can be one file or multiple files. All selected files should be in the same file path location. The following data can be exported: Jobs, Standards, Workspace, Diagnostics and Others (e.g. all files in HunterLab folder, and all pdf reports in Download folder).

- **EXPORT JOB**

This option allows the user to browse and select an existing Job(s) (.ezm) or the current active Job data and copy into a USB flash drive either in **CSV** or **EZM** file format. While exporting into EZM format, the current active Workspace settings can be applied. The color data shown in the Color Data View and the Spectral Data is saved in a CSV file. CSV Files can be exported through **EXPORT > OTHERS**.

Note: In Rev 1.03.0070 and above, the .csv file will be automatically created/updated when a job is saved. CSV files are stored at HUNTERLAB > CSV FOLDER.




Figure 28. Export Current Job

- **EXPORT STANDARD**

This option allows the user to browse and select an existing Standard(s) in the database and copy into the USB flash drive as a file (.std).




Figure 29. Export Standard

- **EXPORT WORKSPACE**

This option allows the user to browse and select an existing Workspace(s) in the database and copy into the USB flash drive as a file (.wsp). To use the above functions, a USB flash drive must be present in the port.

- **EXPORT PDF**

This allows the user to select a PDF file from the Downloads folder to export. Switch to the Download folder in the dropdown list and then select the pdf files to export.



Figure 30. Select Download Folder for PDF File Export

Jobs > Data Management > Email

Saved Jobs can be emailed if there is an active internet connection. When the **EMAIL** option is clicked, the following screen is shown prompting the user to browse and select a user and enter the recipient mail address. You can email any file in HunterLab folder as well as in the downloads folder. Data can be one file or multiple files, e.g. csv file in **HUNTERLAB > CSV FOLDER**, pdf reports in Download folder. In Rev 1.03.0070 and above, the .csv file will be automatically created/updated when a job is saved. CSV files are stored at **HUNTERLAB > CSV FOLDER**. These .csv files can be emailed.




Figure 31. Enter an Address to Email a Job

- **EMAIL SETTINGS**

Click **MAIL SETTINGS** button to configure the SMTP mail server configuration (Port, Server) as shown below. The mail settings configuration is mandatory to enable the mail job feature in the application. When done, press **SEND**.

User's Manual for Aeros and EasyMatch Essentials v 2.0



Figure 32. Enter SMTP Mail Server Information

- **JOBS > DATA MANAGEMENT > DELETE.** The Delete function will allow deletion of Jobs, Standards, Workspace, Diagnostics and others. Data can be one file or multiple files. All selected files must be in the same file path location. In addition, one can delete PDF files from the Downloads folder. .
- **JOBS > DATA MANAGEMENT > BACKUP/RESTORE.** The **BACKUP** function will copy the entire Hunter Lab folder to a thumb drive. **RESTORE** enables the user to copy the backup folder of a thumb drive and upload to the Aeros. .

JOBS > HELP

To access the onboard manual, use Jobs: Help. **NOVICE HELP** can also be enabled under **PREFERENCES > GENERAL**.

JOBS > ABOUT

The **ABOUT** menu provides information about HunterLab and the current software version.




Figure 33. Job > About the Software

To update the software version from a USB flash drive, download the upgrade onto the flash drive. Then insert the USB flash drive into the port on the front of the instrument.



Figure 34. USB port on instrument

Open the **JOBS > ABOUT** menu and press **UPDATE** to continue. After update, open Essentials and it will prompt to enter or create an Administrator Account. If needed, you can edit this account in User Manager later.



User's Manual for Aeros and EasyMatch Essentials v 2.0

For detailed information on firmware, enable a CMR and more, please press the **INFO** button on the screen.




Figure 35. Instrument

Tool Bar: Workspace & System Settings

Workspace Icon



Under the **WORKSPACE & SYSTEMS SETTINGS**, the following tasks can be accomplished:




Figure 36. Workspace Parameters

Workspace > Color Scales

- Color Scales provide four tabs in which the **SCALES, INDICES, DIFFERENCES** and **ILLUMINANT/OBSERVER (ILL/OBS)** can be configured.




Figure 37. Color Measurement Scales

- The Scales Tab shows the five scales available for measurement. Select the absolute scale or color difference scales (if a standard is selected). Press **APPLY** and begin to read your samples.

- The **ILLUMINANT/OBSERVER** tab displays combination selections for these parameters. To see all of the choices, you can scroll through the selections by viewing the screen.




Figure 38. Illuminant/Observer Configuration

- To select indices, go to the **INDICES** tab and check the corresponding box on the right side. Multiple selections are available. To see more choices, the screen can be scrolled. Press **APPLY** to continue. **CUSTOM INDICES** allows user to add reflectance value at a selected wavelength (400-700) as a indice or to customize **HCCI**.



Figure 39. Index Configuration




Figure 40. Custom Indices

User's Manual for Aeros and EasyMatch Essentials v 2.0

- To select differences, go to the **DIFFERENCES** tab and check the corresponding box on the right side. Press **APPLY** to continue.

Color Scales			
Scales	III/Obs	Indices	Differences
dE			<input type="checkbox"/>
dE*			<input type="checkbox"/>
dE CMC			<input type="checkbox"/>
dE* 2000			<input type="checkbox"/>

Defaults Apply Cancel

Figure 41. Color Measurement Differences

Table 3. Overview of Color Measurement Parameters for EZ View, Color Data Table, Trend Plot & Color Plot

Illuminant	Observer	Scales	Differences	Indices	View Options
D65	2/10	CIE Lab	dL*a*b*	457nm Brightness	Pass/Fail ¹
C	2/10	CIE LCh	dL*C*h	Tint E313	Tolerances
F02	2/10	Hunter Lab	dLab	WI E313	Time ³
D50	2/10	XYZ ¹	dXYZ	Y Bright	Date ³
D55	2/10	Yxy ¹	dYxy	YI D1925	Trace Range 1 ²
D75	2/10		dE	YI E313	Trace Range 2 ²
F07	2/10		dE CMC	Z%	Trace Range 3 ²
F11	2/10		dE* 2000	SCAA/G	Trace Range 4 ²
TL84	2/10		dE*	SCAA/C	Auto Range ²
ULT 30	2/10			BCU	Display: Line ²
ULT 35	2/10			HCCI	Display: Point ²
	2/10			Custom ⁴	Zoom
				My, Mc, dM, Tint Strength	Average ²
					Std. Deviation ²
					Meas per Display ²

¹Not Available on Color Plot, ²Trend Plot Only, ³Color Data Table Only, Custom Indices⁴ means selection of a reflectance value at a selected wavelength

Workspace > Read Options

READ OPTIONS > OPTIONS

Shows a dialog box to configure **AVERAGING**, **CONTINUOUS READ INTERVAL**, **AUTO SAVE**, **INDEX BIAS CONFIGURATION**, **SAMPLE NAME**, and **STANDARD NAME**. The Read command performs the operation based on the configured options.




Figure 42. Read Options

Average

Select the number of readings to average to produce the final measurement. The total number of readings to be averaged can be no less than two. Press **APPLY** to close the screen and press **READ** to initiate Readings.




Figure 43. Reading and Averaging

Once the **READ** button is pressed, the instrument will display a unique dialog box to **READ AND AVERAGE** the readings. The second reading is taken using the dialog box button, **READ**. Once all the readings are taken, press **AVERAGE** to obtain the results. Average and Continuous Read are mutually exclusive.

Continuous Read Interval

This feature performs measurements continuously. In **CONTINUOUS READ** mode, measurements are initiated and stopped using the **READ** Button. The minimum value of the Read interval is 3 seconds and it will read as fast as it can update. When in **CONTINUOUS READ** mode, the Read Button is enhanced with a checkmark.

When taking measurements, the Read button is greyed out. When waiting to take the next measurement, the Read button turns green.

Name	L*	a*	b*	Distance (mm)	Height (mm)
Blue Standard	29.36	10.73	-27.85	83.6	5.4
Sample13	29.70	10.59	-27.59	82.8	6.2
Sample12	30.18	10.48	-27.43	82.9	6.2
Sample11	29.48	10.67	-27.74	82.8	6.2
Sample10	29.78	10.58	-27.60	82.8	6.3
Sample9	30.04	10.53	-27.52	83.0	6.1
Sample8	29.37	10.69	-27.78	82.7	6.4
Sample7	29.94	10.55	-27.53	82.8	6.3
Sample6	29.91	10.57	-27.60	83.1	5.9
Sample5	29.39	10.70	-27.79	83.1	6.0
Sample4	30.02	10.52	-27.49	82.8	6.3
Sample3	29.84	10.60	-27.63	83.2	5.9
Sample2	29.67	10.62	-27.66	82.7	6.4
Sample1	29.36	10.73	-27.85	83.6	5.4

Figure 44. Continuous Read

To stop the Continuous Read, press the **READ** button when the button is green.

Auto Save Job

This selection will automatically save a job. Once this feature is selected, a dialog box will be displayed to name the job. If there is no name for a job yet, the file name will default to date, time, instrument and workspace.




Figure 45. Auto Save Job

Index Bias Correction

This option allows the user to input a custom slope and intercept correction for indices. The user can select any Index from the list of applicable indices and input the desired **GAIN** and **BIAS** values. After selecting the required Indices, select the **APPLY** button to save the selected Indices values and update the Views accordingly. The Bias-corrected Indices will be marked with * (e.g.: HCCI *) in the respective view display.

To calculate the slope and bias correction, read a series of samples around the target values of interest. Three methods can be used to provide corrected values:

1. **ONE STANDARD DATA POINT:** In this case, the single data point is compared to the expected value. The Gain remains at 1.0 and the Bias is corrected:

$$\text{BIAS} = \text{EXPECTED VALUE} - \text{MEASURED VALUE}$$

2. **TWO DATA POINTS:** In this case, the two readings are compared to the expected values.

$$\text{BIAS CORRECTION} = \text{EXPECTED VALUE 1} - (\text{MEASURED VALUE 1} * \text{GAIN})$$

$$\text{GAIN CORRECTION} = (\text{EXPECTED VALUE 1} - \text{EXPECTED VALUE 2}) / (\text{MEASURED VALUE 1} - \text{MEASURED VALUE 2})$$

3. **LINEAR REGRESSION:** Create a $y=mx + b$ relationship comparing actual readings to target values, where target values is on the Y-axis and actual readings are on the x-axis. Enter the slope correction under Gain and the intercept correction under Bias.

Select Index Bias Configuration




Figure 46. Slope & Bias Correction



Figure 47. Input Gain & Bias

The indices with bias correction include: 457nm Brightness, BCU, HCCI, SCAA, Tint E313, WI E313, Y Brightness, YI D195, YI E313, Z%, and custom indices.

Prompt For Sample/Standard Name

Select this feature to input the Sample (or Standard) name manually during the measurement cycle so that the Sample measurement will be inserted with the specified name. If this option is not selected, the Samples will be inserted with the specified default sample name suffixed with the auto incremented index number. Press **APPLY** when done.




Figure 48. Prompt for Sample (Standard) Name




Figure 49. Input Sample Name

READ OPTIONS > MEASUREMENT CONFIGURATION

Enables the motion of the turntable, selects the height of the measurement and measurement time.




Figure 50. Measurement Configuration

Select **TURNTABLE MOTION** to enable rotation.

When **AUTO HEIGHT** is selected as a default, the sensor head will adjust its position to read samples automatically for each sample measurement based on the sample height detected by sensor.

FIXED HEIGHT will enable the user to select a height and move sensor head to a fixed location. No matter what sample height is measured by sensor, the sample height that user entered here will be used to do color calculation of each samples. This is usually applied when sensor could not detect sample height very well, e.g. sample has high absorption at 880mm.

The **CREATE SAMPLE PROFILE** button will allow the Aeros to look at the specimen and select the optimum fixed distance from the turntable. Each sample measurement under **CREATE SAMPLE PROFILE** mode will use the height detected by sensor for color calculation. If the detected height is out of preferred sample range (+/-0.5in), it will prompt the error message instead of moving sensor head.

When **CREATE SAMPLE PROFILE** is pressed, the Sensor will move to the Top position.




Figure 51. Creating a Sample Profile

Place sample on the turntable and press **OK**.




Figure 52. Place sample on Turntable




Figure 53. Sample Profile Complete

When profile is complete, press **OK** to continue. The sample height will be shown next to the **CREATE SAMPLE PROFILE** button.

The **MEASUREMENT TIME** is the amount of time used to average the readings together. The longer the time, the more the averaging.

Workspace > Standard And Tolerances

This command can be used to specify the type of **STANDARD AND TOLERANCES**. Standards can be one of four types: **RETRIEVED FROM DATABASE, PHYSICAL, AD HOC** and **NUMERIC**. A standard that is retrieved from Database has been previously stored. A physical standard is one that has been read as a sample and made into a standard. An Ad Hoc (or working) standard is one that is read at the beginning of a job and becomes the standard for a run. In this case, auto tolerances are recommended. A numeric standard is one that has color measurement values but is not present and cannot be read. A subset of this is the Hitch Standard. All types of standards can apply Hitch.




Figure 54. Tolerances Configuration

- Standard is saved with standard name, standard color value and standard tolerances to the database by pressing the button at the bottom of the screen. When there is a standard applied in a job, you must delete it first if you want to change the standard type (Recall. Physical/Adhoc and Numeric). You can click the **CALC AUTO TOLERANCES** here to calculate the tolerances of standards.
- **TOLERANCES** can be entered manually for a selected scale, index and difference.




Figure 55. Enter Tolerances

User's Manual for Aeros and EasyMatch Essentials v 2.0



Figure 56. Index Selection and Tolerances

- Tolerances will be displayed on the measurement screen if enabled under the **VIEW OPTIONS** for the **COLOR DATA** and the **COLOR PLOT** Screens.

The screenshot shows the 'Color Data Table [D65/10]' interface. It displays a table with columns for Name, L*, a*, b*, Distance (mm), Height (mm), Turntable, BCU, dBCU, HCCI, dHCCI, and SCI. The table lists various color standards and coffee samples with their corresponding colorimetric values and tolerances.

Name	L*	a*	b*	Distance (mm)	Height (mm)	Turntable	BCU	dBCU	HCCI	dHCCI	SCI
Standard_201805 03164247	19.47	10.40	13.32	83.04	20.40	On	0.19		6.18		46.
+ Tolerance	0.89	1.04	1.26	0	0	0	1.2	0	8.5	0	0
- Tolerance	0.89	1.05	1.25	0	0	0	1.2	0	8.5	0	0
Coffee58	29.52	4.82	3.32	82.90	20.54	On	1.27	1.08	8.37	2.20	70.
Coffee57	29.51	4.81	3.32	82.34	21.10	On	1.27	1.08	8.37	2.19	70.
Coffee56	29.51	4.82	3.31	82.35	21.09	On	1.27	1.08	8.37	2.19	70.
Coffee5	29.47	4.83	3.32	82.35	21.09	On	1.27	1.07	8.35	2.17	70.
Coffee4	29.34	4.87	3.36	82.67	20.78	On	1.26	1.06	8.30	2.13	69.
Coffee3	27.36	6.72	5.63	81.79	21.65	On	1.06	0.87	8.21	2.04	65.
Coffee2	18.28	9.11	10.95	84.37	19.07	On	0.05	-0.15	5.29	-0.88	42.
Coffee	19.47	10.40	13.32	83.04	20.40	On	0.19	0.00	6.18	0.00	46.

Figure 57. Indices & Tolerances on CDT

- **PASS/FAIL** based on these tolerances can be used on the **EZ VIEW** as well as **COLOR DATA VIEW**.

- HITCH STANDARDIZATION
 - **Hitch Standardization** is a process by which two or more instruments of similar design can be made to read the same color values on a group of specimens. This process is very useful in expanding the communications of color around the world or between vendor and customer.
 - The process of Hitch Standardization (also known as transfer standardization) involves assigning one instrument to be the reference, or master, unit and mathematically adjusting the secondary, or slave, unit(s) to read the “same” values. In this way, two or more instruments can be hitched together. Hitching a secondary unit to a reference instrument requires that a specimen be read on both units and the values compared and adjusted accordingly. This specimen, known as the hitch standard, is first read on the reference instrument and its values recorded as spectral data or colorimetric (tristimulus) data. The hitch standard is then physically moved to the secondary instrument where it is reread and the values from the reference unit are input into the secondary instrument’s processor.
 - Steps for Hitch Standardization:
 1. **READ** a standard.
 2. Go to **WORKSPACE > STANDARD AND TOLERANCES** and select **HITCH**.
 3. The Standard is displayed under **AS READ**. Enter the **TARGET VALUES**.
 4. Check the **APPLY HITCH BOX** and press **OK**.

Hitch Standard

Hitch Method: Tristimulus Hitch

Hitch Type: Additive, Ratio

Colorimetric Conditions

Illuminant/Observer: D65/10

Scale: CIELAB

Readings

	L*	a*	b*
As Read	94.86	-0.92	2.06
Target	<input type="text" value="94"/>	<input type="text" value="1.0"/>	<input type="text" value="2.0"/>

Apply Hitch

OK Cancel

Figure 58. Hitch Standardization

Table 4. Hitch Standardization By Additive or Ratio Application

	Hitch Standard		Calculate Additive	Calculate Ratio	New Read	Apply Hitch Additive	Apply Hitch Ratio
	Target	Measured	Hitch	Hitch		New Read	New Read
	X	80.27	78.29	= +1.98		=*1.025	70.84
Y	81.00	79.21	= + 1.79	=*1.022	72.25	74.04	73.88
Z	50.71	47.76	= +2.95	=*1.061	46.07	49.02	48.91
			Hitch Factor	Hitch Factor			Hitch Calc
			1.98	1.025290586			70.84+1.98 =72.84*1.025
			1.79	1.022598156			72.25+1.79 =72.25*1.022
			2.95	1.061767169			46.07+2.95 =46.07*1.061

AUTOTOLERANCES SETTING

- Using Tab 4, AutoTolerances are calculated for a Color Scale using CMC. The default values of I:C – 2:1 with auto correction factor = 0.75 and commercial factor = 1. However, these ratios can be modified as needed.



Figure 59. AutoTolerance Configuration

Note: If AutoTolerances are selected, the user cannot manually enter tolerances.



Figure 60. Difference Tolerance Configuration

TOLERANCES

- Tolerances can be entered for **SCALES**, **INDICES** and **DIFFERENCES**.



Figure 61. Indices Tolerance Configuration

Workspace > Views

This option can be used to select the views to be presented in the application. Simply check on the box of the screen needed. Press **APPLY** to save one or all of the screens. The default screen is the **COLOR DATA TABLE**. To navigate between screens once the selections have been applied, use the **VIEW FLIPPERS** on the left and right of the screen.




Figure 62. Workspace Views

Workspace > New Workspace

- This allows the user to create a new workspace. A warning is shown to make sure that the current Job is saved. All settings in the previous workspace will be loaded in the new workspace.



Figure 63. New Workspace




Figure 64. Name the New Workspace

System Settings: Standardization

From the **TOOLS** menu select **STANDARDIZATION**. You can also press the Standardization key in the Status bar as a shortcut. The sensor will move to the top to begin.




Figure 65. Begin Standardization

- **READ BLACK GLASS:** Attach the standardization box to the sensor. Then attach the black glass and press **READ** to continue.




Figure 66. Attach Standardization Box with Black Glass

- **READ WHITE TILE:** Remove the black glass and attach the calibrated white tile to the standardization box. Press **READ** to continue.




Figure 67. Change to White Tile

- If the user has selected the **GREEN TILE READING DURING STANDARDIZATION** under **WORKSPACE > PREFERENCES**, the user will be prompted to read the green tile too. The target values are the ones on the back of the green tile.



Figure 68. Green Tile Reading

- Remove the calibration box when standardization is complete. The instrument is ready to read samples.




Figure 69. Standardization is Complete.

System Settings > Diagnostics

Three performance diagnostics and EasyCert are included with software version 1.05.0090 and higher. The three performance diagnostics are REPEATABILITY, GREEN TILE and AUTO DIAGNOSTICS. The EasyCal™ programs under “Validate” offers instrument qualification and performance validation for end-users to self-certify their Vista with traceable standards.



Figure 70. Performance Diagnostics Menu

TESTING THE AEROS FOR COLORIMETRIC REPEATABILITY

The Repeatability Test assesses how consistently the instrument can measure color. To begin the sample pan should be free of samples and obstacles. Click **NEW** button to start the repeatability test and the user is prompted to press **OK** to standardize. All sample readings must be within the tolerances to pass the test.



Figure 71. Set up for Colorimetric Repeatability.




Figure 72. Standardize the Instrument.

- Read the black glass and then the white tile.




Figure 73. Read White Tile

User's Manual for Aeros and EasyMatch Essentials v 2.0




Figure 74. Begin Repeatability Readings on White Tile

- Leave the white tile at the port and press **OK** to initiate the repeatability readings. The white tile is read 30 times and the individual results reported. A table of the difference between the current reading and Standard is shown after every measurement. By comparing each reading to the tolerance, a Pass/Fail assessment is shown.

Result criteria with Tolerance: $dE^* \leq 0.025$								
Name	Result	X	Y	Z	L*	a*	b*	dE*
Standard		81.48	86.62	91.41	94.58	-1.24	1.05	
Sample1	Pass	81.50	86.63	91.43	94.58	-1.24	1.05	0.01
Sample2	Pass	81.51	86.65	91.44	94.59	-1.24	1.05	0.01
Sample3	Pass	81.52	86.66	91.45	94.59	-1.24	1.05	0.02
Sample4	Pass	81.53	86.66	91.46	94.59	-1.24	1.05	0.02
Sample5	Pass	81.53	86.67	91.47	94.60	-1.24	1.05	0.02

Figure 75. Repeatability Readings with Pass/Fail

When all 30 readings have been made, the final test result is shown and saved automatically. To print the results, press the **PRINT** button or click to **OPEN THE FILE** and then **PRINT**.

Name	Result	X	Y	Z	L*
Standard		80.44	85.67	91.85	94.17
Sample1	Pass	80.43	85.67	91.84	94.17
Sample2	Pass	80.43	85.67	91.83	94.17
Sample3	Pass	80.43	85.67	91.83	94.17
Sample4	Pass	80.43	85.67	91.84	94.17
Sample5	Pass	80.43	85.67	91.84	94.17
Sample6	Pass	80.43	85.67	91.84	94.17
Sample7	Pass	80.43	85.67	91.84	94.17

Figure 76. Diagnostics Repeatability Test Results

READING THE GREEN TILE

This test requires entry of the target values for the Green Tile.




Figure 77. Input Target Values for Green Tile

Once the target values have been entered, press **NEXT**. Standardize the instrument and attach the Green Tile. Press **OK** to continue.




Figure 78. Attach the Green Tile

Ten readings are taken and compared to the tolerance as an average. This test is then automatically saved and can be printed by pressing **PRINT**.



Figure 79. Green Tile Readings

RUN AUTO DIAGNOSTICS

Auto Diagnostics is for use by the service department at HunterLab and not recommended for customer use. It runs all tests and detailed readings for short term repeatability and green tile measurements are available by opening the CSV file.

VALIDATE

Aeros Essentials offers instrument validation options for end-users who wish to self-certify their color measurement instrumentation with traceable color standards. Standards are available in individual or three-sample sets, representative of the end user's working color range. Each standard is supplied with a Certificate of Analysis with traceable values and uncertainties.

For more information, please contact HunterLab.

ADVANCED TESTS

Advanced Tests are primarily for use by HunterLab's Service Department. However, this menu supports the upload of Printer Drivers, Parking of the sensor for Shipping, and remote support setting.




Figure 80. Advanced Menu

- **SENSOR HEAD POSITION:** To test the movement of the sensor head up or down. To prepare for shipping, use the **PARK FOR SHIPPING** command under System.
- **READ SIGNAL, DARK, ZERO:** This function will enable the Service Department to determine proper performance of the instrument. These tests can be performed as a single measurement or a continuous loop. The **SIGNAL DATA** for the white tile are shown in the next figure.
- The **TURNTABLE** can be tested next by pressing **START**. Press again to stop the turntable.
- **ENABLE LOG:** Once check enabled, this feature records the instrument actions for tracking purposes. When complete the user returns to this screen and exports the data to a thumb drive. Once the data export is completed the data size becomes '0' again.
- **STANDARDIZE:** Initiates standardization from the Diagnostics screen.
- **MEASURE:** Initiates the measurement of a sample from the Diagnostics screen. The reflectance spectra data will be showed in this dialog.
- **PARK FOR SHIPPING:** Move sensor down to secures the optics for transportation.

User's Manual for Aeros and EasyMatch Essentials v 2.0

- **PRINTER DRIVER:** To upload a new print driver or apk file, download the apk file needed from the internet onto a flash drive. Place the flash drive into the instrument (front port) so that it can access the list of apk files. Select the driver to upload and press **OK**.




Figure 81. Insert USB with Printer Driver



Figure 82. Select Printer Driver




Figure 83. Updating Printer Drivers

- The Aeros will install the new printer driver and it will be then available to use.




Figure 84. Printer Driver Installed




Figure 85. Printer Page

- **RESTART COMM** can be used to reset the ethernet communications for EasyMatch QC.
 - SUPPORT REGION

Note: Your instrument must be connected to the internet.

- a. Select **SUPPORT REGION**. A dialog is displayed with three options for region selection. Select **USA**, **EUROPE** or **ASIA PACIFIC** and the URL string is modified accordingly. Press **OK** to continue.




Figure 86. Select Region

- b. Select **RESTART SUPPORT** to view the Netops Host Screen. From the top right side of the Netops Host screen, select the **3 DOTS**. From the list menu, select **RESTART**.




Figure 87. Remote Access Screen

- c. To ensure that your application is successfully restarted, make sure that you see the message **WEBCONNECT: 'HUNTERLABS' ONLINE**. If this message does not appear, please contact our support team. To exit press the arrow at the bottom of the screen.




Figure 88. WebConnect to HunterLab

PREDICTIVE TESTS

HunterLab Predictive Diagnostic is designed to monitor the software and hardware components of the Agera. Predictive Diagnostic is used to capture different low-level and user-initiated data during normal operation. Following are predictive diagnostics features that have in Aeros 1.05 and above.

1. In **WORKSPACE MENU > DIAGNOSTICS > PREDICTIVE**.
2. Setup reminder interval days for the white (repeatability) and green tiles. Select the number of days for the reminder. Select to **DISABLE/ENABLE** the **TEST EXPIRY ALERT**.
3. Select a trend and press **SHOW** to display the data trend for repeatability, green tile, or monitor channel. Press **APPLY** to view the data trends over time. Export the predictive diagnostics data in csv files to thumb drive.

User's Manual for Aeros and EasyMatch Essentials v 2.0




Figure 89. Predictive Diagnostics

In each trend plot, first select the **TIME RANGE**, then click **SHOW** to display the data. In the plot, click each data point to get the details showing in the right side.




Figure 90. Monitor Channel

4. Select **EXPORT** to send the predictive diagnostics data to a thumb drive. Select the data to export and press **OK**. This data will be sent in CSV format for use in a spreadsheet.




Figure 91. Predictive Test Options

There are three types of predictive diagnostic files.

Diagnostics Data: Records all of the diagnostics tests (i.e. White Tile Repeatability and Green Tile Tests).

Standardization Vector Log: Records raw data from the sample and monitor channel during each standardization.

Initial Data: The original raw data from the sample and monitor channel. This data should not be modified by users.

5. Warning Messages – Collect the following raw data:

Test	Sample to Measure	Warning alert
Standardization vector data		
Sample channel signal Data	Black Glass	Max BOS is above 700.
Monitor Channel Signal Data	White Tile	Max monitor data is below 21000.
Sample Channel Signal Data	White Tile	Max monitor data is below 21000
XYZ Difference	Green Tile	Between 0.25 – 0.3
Service Date	Green Tile	Within 1 month

Once the Disable Error/Warning Alerts is unchecked and applied in Workspace Menu > Diagnostics > Predictive, the info button in the tool bar will list all of the existing warning and error messages. It will be labeled with a different colored dot – Red dot for errors, a yellow dot for warnings and no color for no error or warning,

System Settings > Preferences

PREFERENCES > GENERAL

This allows the user to set preferences to: Load the last used workspace and job, Set standardization time interval, Set screen brightness and date/time, Enable novice tooltip, Enable application security, Use last login credentials, Configure and enable network data export, and Configure network settings.




Figure 92. System Settings> Preferences> General Page

- To **LOAD THE LAST WORKSPACE AT STARTUP** check this box and press **APPLY**.
- To **LOAD THE LAST JOB AT STARTUP**, check this box and press **APPLY**.

User's Manual for Aeros and EasyMatch Essentials v 2.0

- The **STANDARDIZATION TIME INTERVAL** is a useful reminder to restandardize. Press **APPLY** to set the new interval. When the time has lapsed, a prompt to restandardize will be displayed before measurements can be taken.
- Set the screen **BRIGHTNESS** using a sliding scale and press **APPLY**.
- Set the **DATE** and **TIME**, time zone, and format use the **ADJUST CLOCK** feature.
- **LANGUAGE** can be selected here along with a keyboard change. Languages supported include English, Chinese, Japanese and German.




Figure 93. Language & Keyboard Selection

- **ENABLE NOVICE TOOLTIPS** by checking on the box. Once enabled, on screen tips are displayed for 3 seconds. To display again, roll over the **LIGHTBULB ICON** on the lower right part of the screen.




Figure 94. Enable Novice Tool Tips



Figure 95. Example of Novice Tool Tip

- **ENABLE APPLICATION SECURITY.** This selection is available after the User Manager has been set up. Please refer to the **SYSTEM SETTINGS > USER MANAGER** for more information.
 - When this is selected, the application will require valid login credentials at startup. On successful login, the user name will be shown in the status bar. If **USE LAST LOGIN CREDENTIALS** is checked, the user will be automatically logged in on subsequent startups.
- **GREEN TILE CHECK** adds the green tile reading to the standardization process.
- **AUTO DIM** darkens the screen after Essentials idle 15min. To disable, check the box **DISABLE AUTO DIM**.
- To **CONFIGURE AND ENABLE THE NETWORK DATA EXPORT** and **NETWORK SETTINGS**, please connect the Aeros to a computer as described in **SPECIAL FUNCTIONS** (Chapter 7). In Network settings, one can use an Ethernet cable.

PREFERENCES > PRINT

This page allows the user to configure:

- The **READINGS** to print
- The option to **PREVIEW** before print.
- **PRINT REPORT TITLE** and **LOGO**. To apply a logo, import the logo first to the HunterLab folder.
- Orientation of the report (**Portrait or Landscape**) orientation.
- To save changes, press **APPLY**.




Figure 96. System Settings> Preferences> Configure Print Page

System Settings > User Manager

Security can be enabled on the Aeros to ensure that operators cannot modify/delete folders or files and limit their functionality. An administrator is identified to set up the users/groups with selected privileges.

- To begin, go to **WORKSPACE > USER MANAGER** to set up **CREATE ADMINISTRATIVE GROUPS** first followed by **CREATE USER GROUPS**.




Figure 97. Create a Group.

- Once the groups have been established, then individual users with **USER NAMES** and **PASSWORDS** can be setup for both Administrator and User Groups.



Figure 98. Setup Administrative & General Users

User's Manual for Aeros and EasyMatch Essentials v 2.0

- Users in Administrative Groups have all features enabled. For Users in User Groups, Privileges can be setup as shown below. Press **UPDATE PROFILE** to complete.




Figure 99. User Privileges

- If a printer is attached, you can **PRINT** a list of Privileges selected.
- To complete enabling security, go to **WORKSPACE > PREFERENCES** and enable security on the right side.




Figure 100. Enabling Security

- After enabling security, each user must enter a name and password when logging into the Aeros. For convenience, the user can check the box under **WORKSPACE > PREFERENCES > GENERAL** to use the **LAST LOGIN CREDENTIALS**.



Figure 101. Login Credentials

User's Manual for Aeros and EasyMatch Essentials v 2.0




Figure 102. Enable Last Login Credentials

- If needed, the administrative user can delete groups / users and reset passwords of all Groups & Users.
- After **ENABLE APPLICATION SECURITY** is checked, the Logoff feature will be listed in Job menu. User can click **LOGOFF** to exit Essentials.




Figure 103. Log Off

Toolbar: Options (Views)

Views Icon



VIEWS are selected using a dialog box under **WORKSPACE**. Simply check on the box of the screen needed. Press **APPLY** to save one or all screens. The default screen is the Color Data Table. To navigate between screens once the selections have been applied, use the **VIEW FLIPPERS** on the left and right of the screen.



Figure 104. Workspace > select Views.

Once the views have been selected, then **VIEW OPTIONS** are available to add additional information to the screen. Each screen has a unique set of options associated with it.

Views > EZ View

This view provides a simple display of **STANDARD** vs. **SAMPLE** and **PASS/FAIL** results.

Name	Standard_20...	Sample
L*	96.54	----
a*	-0.73	----
b*	5.63	----
dE*		----
WI E313 [C/2]	66.15	----

Figure 105. EZ View Display

- **VIEW OPTIONS** includes the selection of **NO COLOR SCALE**, **PASS/FAIL**, **PRECISION** and **FONT SIZE**.




Figure 106. EZ View Options

Views > Color Data Table

The **COLOR DATA TABLE** view shows **COLOR SCALE**, **COLOR DIFFERENCE**, and **INDEX** data for the **STANDARDS** and **SAMPLES** in the job.

The screenshot shows the "Color Data Table [D65/10]" interface. It features a table with the following data:

Name	L*	a*	b*	dE*	WI E313 [C/2]
Standard_202009 23163447	96.54	-0.73	5.63		66.15

The interface also includes a top toolbar with icons for refresh, up, eye, settings, and menu. A status bar at the bottom shows "Standardized", "Job: 20201123_115903_ARIS00006_Default*", and "Workspace: Default".

Figure 107. Color Data Display

- Options such as **TOLERANCES**, **DATA ORDER**, **INTERVAL**, **DATE**, **DISTANCE**, **TURNTABLE MOVEMENT**, **HEIGHT**, **TIME**, **USER NAME**, **SENSOR NUMBER** and **PASS/FAIL** can be selected for viewing using the **VIEW OPTIONS**..



Figure 108. Color Data Screen: View Options Needs new screen.

- A long press on the **SAMPLE** can enable the user to make the sample into a Standard, change the name or delete the reading.




Figure 109. Changing a Sample into a Standard

- To delete a Sample (or Standard), select **DELETE** and then confirm the action.




Figure 110. Delete the Sample Measurement

A long press on the **STANDARD** will enable the user to edit or delete the Standard. Edit opens the Standard and Tolerances dialog box. Delete will delete the standard from the current workspace.



Figure 111. Edit Standard Name or Delete the Standard

Views > Spectral Data Table

The **SPECTRAL DATA TABLE** displays percent reflectance or absorbance values for each selected measurement at the wavelengths being measured.



Figure 112. Spectral Data Table

- Selections include **ABSOLUTE OR DIFFERENCE, REFLECTANCE OR K/S**. Enter the start and stop wavelength, the interval and the precision and press **OK** to continue.




Figure 113. Spectral Data Table Options

Views > Spectral Plot

- This view provides a plot of wavelength vs. spectral measurement parameter.



Figure 114. Spectral Plot View

- Press **CLEAR ALL** to remove all the samples to display. Press **SELECT ALL** to enable display of all samples. To select an individual sample, click on the respective Sample in the list located on the right edge of the screen.
- The Samples List is paginated. Click the **LEFT** and **RIGHT ARROW** buttons below the samples list to navigate between pages.
- Press and hold on the left/right page number arrows under the sample list to show a small dialog box. This dialog allows you to select the number of records per page to display and the page number to display.




Figure 115. Spectral Plot Options

- **Spectral Plot Options:** There are three choices for spectral plot options:
 - **K/S** – mathematical calculation based on reflectance and determined at each wavelength for the standard and sample.
 - **REFLECTANCE** - Displays the reflectance value at each wavelength.
 - **% STRENGTH** – Percentage of the ratio of the K/S of the sample to the K/S of the standard.
- Uncheck the **OPTIONS > SHOW BACKGROUND**, to display the plot with white background color.

- Check **OPTIONS > AUTO RANGE** to automatically scale the contents to fit. If **AUTO RANGE** is not selected, then enter the **Y-** and **X-axis** range to display.

Views > Trend Plot

This tool can be used to study the trends in production and identify color variations. There are four parameters of color measurement (three scale values and optional indice) which can be represented in four traces. If a sample point is selected in one of the four traces, it is highlighted in blue. The name is shown at the bottom right hand corner of the View. The **AVERAGE** and **STANDARD DEVIATION** can be shown as per the view configuration settings.




Figure 116. Trend Plot

- **VIEW OPTIONS** for the **TREND PLOT** includes the type of display, the statistics and the number of readings per display.




Figure 117. Trend Plot Options

- **VIEW OPTIONS > TRACES** set the ranges for the traces or allow selection of **AUTO RANGE**. Trace 1 to 3 uses the current Color Measurement Scale and Trace 4 will allow for measurement of differences or an index. The user can select which Traces to view and set control limits as a percent.




Figure 118. Trend Plot Traces

Views > Color Plot

This shows the sample location in two-dimensional Color Space with respect to the standard for difference measurements or the samples in absolute measurement. For differences, the standard is the center point of the plot and the samples are plotted separately on the graph.

The displayed samples are shown in a list box on the right of the screen. The color plot can be zoomed, and the data points can be viewed in detail.

Press and hold on the left/right page number arrows to show a small dialog box. This dialog allows you to select the number of records per page to display and the default page number to display.

Figure 119. Color Plot View



Figure 120. Color Plot Options

The tolerance plot is available in rectangular and elliptical color space. The **PASS/FAIL** sample points are shown in green and red when in difference mode, respectively. In Absolute Mode, they are shown in green.

Recall Standard

Recall Icon



The Tool Bar offers a quick way to recall standards and samples using this icon.

Click **RECALL** option in Data Management, a dialog will be displayed where the user can recall the measurements from the database by **SELECTING THE TYPE**:

- Show The Standard Associated With A Specified **CATEGORY**.
When this options is selected, the standards list is filled with the standards associated with the selected Category from the **SELECT STANDARD CATEGORY** list. All matching samples are filled into the sample list. To narrow the list, the user can select the samples only associated with the selected standard.
- Show All The Measurements In The Selected **JOB**
When this option is selected, the standard and samples for that job are listed.
- Show The **STANDARDS/SAMPLES** In The **CURRENT WORKSPACE**
When this option is selected, the list is filled with the standards and samples associated with the current workspace. To narrow the sample list, the user can select the samples only associated with the selected standard.
- Show The **IMPORTED STANDARDS**
When this option is selected, the standards list is only filled with the standards imported into the database. After selecting the measurements click **RECALL** placing them into the active Job.




Figure 121. Recall Measurements

User's Manual for Aeros and EasyMatch Essentials v 2.0




Figure 122. Select Standard by Type

Select the Type and then press Recall to continue.

Special Functions

Auto-Exporting Data Through A Network Connection

Connect Aeros to a Network. You can connect Aeros to a network hub using the Ethernet cable. The computer must be connected to the same network as the Aeros.

CONNECT TO A NETWORK HUB USING AN ETHERNET CABLE

- Hardware needed: Ethernet cable plugged into the back of the Aeros and the other end to a network hub.



Figure 123. Ethernet Cable

- To connect Aeros to network, go to **WORKSPACES > PREFERENCES** and Select **CONFIG NETWORK SETTINGS**.




Figure 124. Preferences (General) > Network Settings

User's Manual for Aeros and EasyMatch Essentials v 2.0




Figure 125. Configure Ethernet

- Select **CONFIGURE ETHERNET SETTINGS**.
- Check **USE DHCP FOR ETHERNET CONFIG**. Please write down the IP address showing in the Ethernet Setting dialog. You can also check the IP address of Aeros in **JOBS > ABOUT > INFO**.




Figure 126. Select DHCP

- Go to **WORKSPACE > PREFERENCES** and select **AUTO NETWORK DATA EXPORT MEASUREMENT** using a check and select **CONFIG**. Choose **ACT AS SERVER** and **PORT** number as **11111**. You can also choose a delimiter to mark your data.




Figure 127. Auto Export Measurement




Figure 128. Network Data Export

- Configure the computer with the following settings:
 - ❖ Set computer as **CLIENT**.
 - ❖ Enter the **IP ADDRESS** of Aeros as recorded above.
 - ❖ Set the port number as **11111**.
- After all have been set, the data is ready to be exported from Aeros to the computer.




Figure 129. Data Export

Auto-Exporting Data Via Direct Connection Between Aeros And A Computer

Ethernet cable is plugged into the back of the Aeros and the other end is connected to the computer. Ethernet adapter USB can be applied here if the computer does not have available Ethernet port.

- Materials Needed: Ethernet cable and Ethernet adapter to USB can be applied here if the computer does not have an available Ethernet port. Hardware needed: Ethernet cable and Ethernet adapter to USB can be applied here if the computer does not have available Ethernet port.



Figure 130. Ethernet Cable



Figure 131. Ethernet to USB Adapter

- Connect Aeros to Computer:
 - Plug Ethernet cable into RJ-45 Ethernet connection at rear of Aeros.



Ethernet Port

Figure 132. Rear View of Aeros

- Plug the other end into the Computer or into the Ethernet adapter
- Open Command Prompt in the PC.

Type in **IPCONFIG**, find the right ethernet (in this case, it is **Ethernet adapter Ethernet**) and write down **AUTOCONFIGURATION IPV4 ADDRESS** as well as the **SUBNET MASK**.




Figure 133. Command Prompt ipconfig

- Configure the Aeros

Open Aeros Essentials, go to **WORKSPACES > PREFERENCES > CONFIGURE NETWORK SETTINGS**. First, select the Ethernet configuration. Uncheck **USE DHCP FOR ETHERNET CONFIG**. Type in **IP ADDRESS** and **SUBNET MASK** manually. The IP address here should be same as the autoconfiguration IPv4 Address in the PC, except changing the last number. The Subnet Mask is the exact same. Restart Aeros to get network setting applied.

User's Manual for Aeros and EasyMatch Essentials v 2.0




Figure 134. Configuration Parameters for Ethernet

- Press **APPLY** on the Ethernet Configuration and then **APPLY** on the Preferences Page to complete.
- Turn the instrument **OFF** and then back **ON**.
- Go to **PREFERENCES** and select **AUTO NETWORK DATA EXPORT**.




Figure 135. Read Options > Auto Export Measurements

- For a direct connection between Aeros and data collection computer, set up the Aeros as a **CLIENT**.
- Enter the **COMPUTER IP ADDRESS** here, in this case **169.254.113.144** and the **PORT** as **10001**. Press **APPLY** on the screen to continue.




Figure 136. Read Options Export

- Aeros is now ready to send data.
- Configure the Computer:
 - Connection configurations differ depending on data collection software. The data collection computer will be set up as a Server.
 - Connect as follows:
 - ❖ Set computer as **SERVER**.
 - ❖ Enter the computer IP address **169.254.113.144**
 - ❖ Put the port number as **10001**
- Send Data from the Aeros:
 - Configure the **COLOR DATA TABLE** with the color scale and parameters to be measured.




Figure 137. Data Output

Tips & Tricks: Recover Unsaved Measurement Data

- In the case where the application is closed unexpectedly, the data is temporarily stored in a table along with the Job details. When the application restarts, a prompt allows the user to recover the data.




Figure 138. Recover Data

- If the user answers **YES**, all measurements are recovered into a new job or appended to a saved job.

HunterLab File Service Package

The HunterLab File Service is a customized background service which provides the network storage facility for Essentials-AEROS to back-up a File or Folder to a networked PC. This package contains:

- HunterLab File Service Installer (FileServiceInstaller.exe)
- A package file **HLFSPACKAGE.PKG**.

TO INSTALL THE FILE SERVICE PACKAGE

- Copy the above installation package files into a networked PC.
- Run the executable **FILESERVICEINSTALLER.EXE** and follow the guided steps to complete the installation process.
- After installation a shortcut for **HUNTERLAB FILE SERVICE CONFIG TOOL** will be created on Desktop. Double click on the shortcut **HUNTERLAB FILE SERVICE CONFIG TOOL**.




Figure 139. File Service Tool

- Select the **ROOT FOLDER** by clicking on the **BROWSE** button. The Essentials Backup will be stored in the configured Root Folder path. Enter the **PORT NUMBER** for the network File service. Click on the **APPLY** button. The File service will be restarted with the new settings.




Figure 140. Select Root Folder

Note: Make sure that the configured port number is added to Exceptions in the firewall. The configured port number in the server must be same at the client side (Essentials-AEROS).

SETTING UP FILE STORAGE FROM AEROS (CLIENT) SIDE

In Aeros Essentials, navigate to **JOBS > DATA MANAGEMENT > BACKUP (or RESTORE)**. The **SELECT ACTION** dialog will be displayed. The user can choose between **USB DRIVE** or **NETWORK STORAGE**. When **USB DRIVE** option is selected, the Backup and Restore operations will be performed into the USB flash drive plugged into the system.




Figure 141. Select USB Option

When **NETWORK STORAGE** is selected, the Backup and Restore operations are performed into a network folder of the specified system where the HunterLab File Service is installed. Click on **NETWORK STORAGE SETTINGS** Button.




Figure 142. Network Storage Settings

In the next screen enter the **IP ADDRESS** and **PORT NUMBER**. Click on **TEST CONNECTION** button to verify the connectivity. Click **APPLY** to save the settings. The saved network settings will be used for the Network Backup and Restore operations.




Figure 143. Network Storage Settings

After successful configuration of network settings, click **BACKUP** (or **RESTORE**) to perform the complete backup of **HUNTERLAB** folder in Essentials-Aeros to the specified network server's folder.




Figure 144. Select Restore

Select the files to be restored.




Figure 145. Select Files to be Restored.

User's Manual for Aeros and EasyMatch Essentials v 2.0




Figure 146. Files Restored

Specifications

The specifications and characteristics of the instrument are given in this chapter. For best performance, the instrument should be placed where there is ample workspace with medium or subdued illumination and no drafts. The operating conditions (temperature and humidity ranges) are given in the Operating Conditions section below.

Note: Do not leave Aeros in an area where temperature or humidity extremes are possible.

Operating Conditions

Storage Temperature (3weeks)	-20°C to 65°C (-5°F to 150°F)
Operating Temperature	4°C to 38°C (40°F to 100°F)
Noncondensing Humidity	10% to 90% for storage 2-% to 85% for operating
Standard Accessories	Calibrated instrument White Tile, Certificate of Traceability, Black Glass Standard, Green Diagnostic Tile, Standards Box, 30.5 cm (12-in) and 15.2 cm (6 in) sample pans, Power Supply, Quick Start Guide, Aeros User's Guide on CD

Physical Characteristics

Weight	23.0 kg (50.0 lbs.)
Dimensions (Height x width x depth)	56 cm x 38 cm x 51 cm (22.0 in x 15 in x 20 in)
Display	Touch Screen, high resolution color, 7in, 1280x800
Maximum Sample Height	65 – 140 mm (2.5 – 5.5 in)
Communications I/O: USB	Connectivity to printer, keyboard, mouse and other peripherals. Front Panel USB: 2.0 bidirectional data export/import via thumb drive.
Ethernet RJ45	Print directly to standalone or network printers; email directly from the instrument; stream data to LIMS and SPC systems.
Remote Access Support	Enabled via internet-based support tool
System Power	100 – 240 VAC, 47 – 63 Hz to universal power supply @ 24 VDC/3.75A 90W

Conditions Of Illumination And Viewing

Light Source	Full spectrum balanced LED system array; 5 year typical LED life
Dual Beam Spectrophotometer	256 element diode array and high resolution, concave holographic grating; Sealed optics
Measurement Principle	Dual-beam Non-Contact Reflectance Spectrophotometer
Measurement Method	Port down, Npn-contact, Rotating platter @ 12 RPM
Area Measured per Rotation	177.25 cm ² (27.5 in ²)


Instrument Performance

Spectral Data	Range: 400-700 nm; Reporting Interval (nm): 10 nm
Spectral Resolution	<3 nm
Spectral Component	Excluded
Effective Bandwidth	10 nm equivalent triangular
Sampling Rate	Continuous at 7 measurements per second
Photometric Range	0-150%
Measurement Duration	User Selectable 5(, 10, 15, 30 seconds)
Measurements per Rotation	35
Maximum Sample Height	140mm (5.5in)
Inter-instrument Agreement	$\Delta E^*L^*a^*b^* < 0.30$ (Avg) on CCSII BCRA Tile Set
Colorimetric Repeatability	$\Delta E^* < 0.025$ Max on White Tile

Measurement

Data Views	Color Data, Spectral Plot, EZ View, Color Plot, Trend Plot, Spectral Data
Illuminants	A, C, D50, D55, D65, D75, F02, F07, F11, TL84, ULT30, ULT35
Observers	2° and 10°
Color Scales	CIE L*a*b*, Hunter Lab, CIE L*C*h, CIE Yxy, CIE XYZ and differences
Color Difference Indices	ΔE^* , ΔC^* , ΔE , ΔE CMC, ΔE 2000
Indices and Metrics	E313 Yellowness, E313 Whiteness, YI D1925, Y Brightness, Z%, 457nm Brightness, Baking Contrast Units, Tint, HCCI, SCAA
Data Storage	>1 million Records; 8GB
Languages	English, Japanese, simplified Chinese
External Software	HunterLab EasyMatch QC, EasyMatch QC-Electronic Records

Regulatory Notice



HunterLab
ISO 9001 Certified

Declaration of Conformity

EU / EMC Directive:	2014/30/EU
Standard to which Conformity is Declared:	IEC 61326-1: 2012 / EN: 2013
Manufacturer:	Hunter Associates Laboratory, Inc. 11491 Sunset Hills Rd, Reston, VA, USA
European Representative: Representative's Address:	Christian Jansen Griesbraeustrasse 11, 82418 Murnau, Germany
Type of Equipment:	Reflectance Spectrophotometer
Model No.:	Aeros

I, the undersigned, hereby declare that the equipment specified above conforms to the Directive(s) and Standard(s) above

Place: <u>Reston, VA, USA</u>	Signature: <u></u>
Date: <u>December 14, 2017</u>	Full Name: <u>Tim Barrett</u>
	Position: <u>Systems Engineer</u>

Features, Accessories & Maintenance

Aeros Maintenance & Safety

The Aeros is engineered to be virtually maintenance free. This section outlines the few parts of the sensor that are to be maintained for the instrument to function properly.

CLEANING THE AEROS

The Aeros is NOT waterproof, but the exterior of the case may be wiped with a damp cloth.

CLEANING THE WHITE TILE

The White Standard is an optical coating and should be handled in much the same way as other optical surfaces. Although the material is very durable, care should be taken to prevent contaminants such as finger oils from contacting the material's surface. If the surface appears lightly soiled, it may be air brushed with a jet of clean dry air. For heavier soil, the material can be cleaned by scrubbing with a soft brush under running water. Blow dry with clean air or allow the material to air dry. If the material is heavily stained, soak with either an extremely mild mix of soap and water, 5% white distilled vinegar, or hydrogen peroxide. Then run under water while scrubbing with a soft brush. Always keep tiles in the Standards box when not in use

CLEANING THE BLACK GLASS AND GREEN TILE

The **Green tile and Black Glass** can be cleaned using a soft nylon-bristle brush, warm water, and laboratory grade detergent such as SPARKLEEN. Wipe the tiles dry using a clean, non-optically brightened, lint free paper towel, or use warm water as a rinse and let stand to air-dry in a couple of minutes.

Note: SPARKLEEN is manufactured by Fisher Scientific Co., Pittsburgh, PA 15219, and may be ordered from them using catalog number 4-320-4. Add one tablespoon of SPARKLEEN to a gallon of water.

The above procedure is particularly useful if the lab area is not clean. If, however, the lab is clean, an equally effective method for occasional tile cleaning is to use IPQ (isopropyl alcohol) sprayed onto a clean, non-optically brightened, lint free paper towel such as a Kim wipe. Wipe tile thoroughly watching for fingerprints and let air dry.

Keep the **Black Glass** in the standards case when not in use to prevent it from becoming scratched or collecting dust. Before standardizing the instrument, check the black tile for scratches and dust. Significant scratches that result in a hazy appearance to the finish may cause standardization to be in error. If the black tile is scratched, call the HunterLab Order Processing Department or contact your local HunterLab representative to order a replacement.

POWER REQUIRED

- Voltage: 100-240 VAC, 3.75A, 47/63 Hz
- Single Phase
- 60 VA maximum
- Installation Category (Over Voltage): II

SAFETY

- Do not view the instrument LED's directly as it may be damaging to the eyes.
- Do not submerge the instrument in water.
- Do not take the instrument apart as there are 'no user serviceable parts' in the instrument.
- Do not disassemble the instrument and attempt to clean the optical components.
- Do not open the instrument or remove any covers except using the instructions given in this User's Manual or under the direction of HunterLab Technical Support.

When You Need Assistance

If you need for technical or sales assistance on applications, troubleshooting, , service, warranty, accessory pricing and more, please contact the office nearest you:

- For the Americas, Support@hunterlab.com
- For Asia, AsiaSupport@hunterlab.com
- For Europe, EuropeSupport@hunterlab.com
- For India, Middle East and Africa, IMEASupport@hunterlab.com
- For all other regions, Support@hunterlab.com

Additionally, our global support website offers 24/7 assistance with a library of information on various color measurement and appearance topics such as applications, instrument operation, and troubleshooting. The HunterLab global support website is located at support.hunterlab.com.

For personalized assistance, go to support.hunterlab.com and locate the [Create A Ticket](#) button on the menu. A subsequent form gathers information on your request for response from our Customer Experience Teams around the globe.

Index

Auto Save Job, 38
Average, 37
Backup/Restore, 29
Color Data screen, 15
Color Data Table, 11, 36, 46, 65, 66
Color Plot View, 71
Color Scale, 15, 66
Color Scales, 33
Continuous Read, 37
Data Management, 23
Delete, 29
Email, 28
Ethernet Port, 8
Export, 27
EZ View, 36, 43, 65
Import, 25
Index Bias Correction, 38
Instrument Power, 8
Jobs, 19, 21
Keyboard & Mouse, 8
Mail Settings, 28
New Workspace, 47
Operating Conditions, 85
Read Options, 37
Recall, 24
Sample Name, 15, 39
Sensor Status, 18
Spectral Data Table, 68
Spectral Plot View, 69
Standard and Tolerances, 42
Standardization, 11
Standardize, 47
Standards & Tolerances, 15
Taking a Measurement, 11
Tolerances, 66
Tools, 47
Trend Plot, 36, 70
USB Connectors, 8
User Manager, 62
View Flippers, 17, 46, 65
View Options, 19, 65, 66, 70
Views, 46
Workspace, 11, 15, 18, 33

Table of Figures

Figure 1. Calibration Box.....	7
Figure 2. Rear View with Connectors	7
Figure 3. USB Port on Rear of Instrument	8
Figure 4. Measurement Screen	11
Figure 5. Sensor Moves to Top	12
Figure 6. Read the Black Glass for Bottom-of-Scale	12
Figure 7. Read the White Tile for Top-of-Scale.....	12
Figure 8. Completed Standardization	13
Figure 9. Standardization Status Bar	13
Figure 10. New Workspace	13
Figure 11. Read Options: Measurement Configuration.....	14
Figure 12. Configure Sample Name	15
Figure 13. Changing, Renaming or Deleting a Sample	16
Figure 14. Edit/Delete a Standard	16
Figure 15. User Interface Screen for Aeros & Essentials	17
Figure 16. Sensor Serial Number	18
Figure 17. Job Menu	21
Figure 18. Open A Job.....	22
Figure 19. Save A Job	22
Figure 20. Workspace & System Settings > Preferences > Print.....	22
Figure 21. Save as PDF.	23
Figure 22. Save PDF to Downloads File.....	23
Figure 23. Download File Contents.....	23
Figure 24. Data Management Menu.....	24
Figure 25. Recall Measurements	25
Figure 26. Import Job.....	26
Figure 27. Import Workspace	26
Figure 28. Export Current Job	27
Figure 29. Export Standard	27
Figure 30. Select Download Folder for PDF File Export	28
Figure 31. Enter an Address to Email a Job	28
Figure 32. Enter SMTP Mail Server Information	29
Figure 33. Job > About the Software	29
Figure 34. USB port on instrument	30
Figure 35. Instrument	31
Figure 36. Workspace Parameters.....	33
Figure 37. Color Measurement Scales	33
Figure 38. Illuminant/Observer Configuration	34
Figure 39. Index Configuration	34
Figure 40. Custom Indices.....	34
Figure 41. Color Measurement Differences	35
Figure 42. Read Options	37
Figure 43. Reading and Averaging	37
Figure 44. Continuous Read.....	38
Figure 45. Auto Save Job.....	38
Figure 46. Slope & Bias Correction	39
Figure 47. Input Gain & Bias	39
Figure 48. Prompt for Sample (Standard) Name	40
Figure 49. Input Sample Name	40
Figure 50. Measurement Configuration	40

User's Manual for Aeros and EasyMatch Essentials v 2.0

Figure 51. Creating a Sample Profile.....	41
Figure 52. Place sample on Turntable	41
Figure 53. Sample Profile Complete	41
Figure 54. Tolerances Configuration.....	42
Figure 55. Enter Tolerances	42
Figure 56. Index Selection and Tolerances	43
Figure 57. Indices & Tolerances on CDT	43
Figure 58. Hitch Standardization	44
Figure 59. AutoTolerance Configuration	45
Figure 60. Difference Tolerance Configuration	45
Figure 61. Indices Tolerance Configuration	46
Figure 62. Workspace Views.....	46
Figure 63. New Workspace.....	47
Figure 64. Name the New Workspace	47
Figure 65. Begin Standardization.....	47
Figure 66. Attach Standardization Box with Black Glass.....	48
Figure 67. Change to White Tile	48
Figure 68. Green Tile Reading.....	48
Figure 69. Standardization is Complete	49
Figure 70. Performance Diagnostics Menu.....	49
Figure 71. Set up for Colorimetric Repeatability.	50
Figure 72. Standardize the Instrument.....	50
Figure 73. Read White Tile.....	50
Figure 74. Begin Repeatability Readings on White Tile	51
Figure 75. Repeatability Readings with Pass/Fail	51
Figure 76. Diagnostics Repeatability Test Results	51
Figure 77. Input Target Values for Green Tile.....	52
Figure 78. Attach the Green Tile.....	52
Figure 79. Green Tile Readings	52
Figure 80. Advanced Menu	53
Figure 81. Insert USB with Printer Driver	54
Figure 82. Select Printer Driver.....	54
Figure 83. Updating Printer Drivers	55
Figure 84. Printer Driver Installed.....	55
Figure 85. Printer Page	56
Figure 86. Select Region	56
Figure 87. Remote Access Screen	57
Figure 88. WebConnect to HunterLab	57
Figure 89. Predictive Diagnostics.....	58
Figure 90. Monitor Channel.....	58
Figure 91. Predictive Test Options.....	58
Figure 92. System Settings> Preferences> General Page	59
Figure 93. Language & Keyboard Selection	60
Figure 94. Enable Novice Tool Tips	60
Figure 95. Example of Novice Tool Tip.....	60
Figure 97. System Settings> Preferences> Configure Print Page	61
Figure 97. Create a Group.....	62
Figure 98. Setup Administrative & General Users	62
Figure 99. User Privileges	63
Figure 100. Enabling Security	63
Figure 101. Login Credentials	63
Figure 102. Enable Last Login Credentials	64
Figure 103. Log Off.....	64

User's Manual for Aeros and EasyMatch Essentials v 2.0

Figure 104. Workspace > select Views.	65
Figure 105. EZ View Display.....	65
Figure 106. EZ View Options.....	66
Figure 107. Color Data Display.....	66
Figure 108. Color Data Screen: View Options Needs new screen.	67
Figure 109. Changing a Sample into a Standard.....	67
Figure 110. Delete the Sample Measurement.....	67
Figure 111. Edit Standard Name or Delete the Standard.....	68
Figure 112. Spectral Data Table.....	68
Figure 113. Spectral Data Table Options.....	68
Figure 114. Spectral Plot View.....	69
Figure 115. Spectral Plot Options.....	69
Figure 116. Trend Plot.....	70
Figure 117. Trend Plot Options.....	70
Figure 118. Trend Plot Traces.....	71
Figure 119. Color Plot View.....	71
Figure 120. Color Plot Options.....	71
Figure 121. Recall Measurements.....	73
Figure 122. Select Standard by Type.....	74
Figure 123. Ethernet Cable.....	75
Figure 124. Preferences (General) > Network Settings.....	75
Figure 125. Configure Ethernet.....	76
Figure 126. Select DHCP.....	76
Figure 127. Auto Export Measurement.....	76
Figure 128. Network Data Export.....	77
Figure 129. Data Export.....	77
Figure 130. Ethernet Cable.....	77
Figure 131. Ethernet to USB Adapter.....	78
Figure 132. Rear View of Aeros.....	78
Figure 133. Command Prompt ipconfig.....	78
Figure 134. Configuration Parameters for Ethernet.....	79
Figure 135. Read Options > Auto Export Measurements.....	79
Figure 136. Read Options Export.....	80
Figure 137. Data Output.....	80
Figure 138. Recover Data.....	81
Figure 139. File Service Tool.....	81
Figure 140. Select Root Folder.....	81
Figure 141. Select USB Option.....	82
Figure 142. Network Storage Settings.....	82
Figure 143. Network Storage Settings.....	83
Figure 144. Select Restore.....	83
Figure 145. Select Files to be Restored.....	83
Figure 146. Files Restored.....	84