


# Accessories

## Accessories

# ColorFlex® EZ



	<b>Port Inserts</b>
04-6623-10	32 mm (1.25 in)
02-6624-05	32 mm (1.25 in) with glass
D02-1010-618	32 mm (1.25 in) with 420nm UV Filter Assembly
04-6623-09	25 mm (1 in)
02-6624-01	25 mm (1 in) with glass
04-6623-08	19 mm (0.75 in)
04-6623-05	13 mm (0.50 in)



04-6622-00	<b>Port Insert for Sample Cup</b> 64 mm (2.5 in) The use of the sample cup enables the measurement of powders, granules, pellets, pastes and liquids. The port insert is designed to hold the sample cup in proper position at the measuring port.
------------	---



04-7209-00	<b>Glass Sample Cup: 64 mm (2.5 in)</b> Optically clear glass cup for sample presentation of liquids, powders, granules, and pellets.
------------	--



02-4579-00	<b>Ring and Disk Set</b> For use in the glass sample cup, used for transparent and translucent liquids on semi-solids where path length must be fixed. Does not include sample cup.
------------	--



04-4000-00	<b>Sample Cup Opaque Cover</b> Provides a light trap to exclude the interference of external light on the sample cup.
------------	--



D02-1010-327

**Footswitch Assembly**

The footswitch is used to initiate sample measurement without the use of hands.

D02-1010-459

**Sample Clamp / Port Forward Stand**

Securely holds the instrument in a port forward orientation. The sample clamp permits to be held flat against the sample port.



L02-1014-582

**Hard-Sided Carrying Case, CFEZ**

Hard-sided Carrying Case for Standard CFEZ

L02-1014-607

**Hard-Sided Carrying Case, CFEZ Citrus**

Hard-sided Carrying Case for Citrus CFEZ

L02-1014-608

**Hard-Sided Carrying Case, CFEZ Tomato**

Hard-sided Carrying Case for Standard CFEZ



A13-1014-294

**USB Flexible Keyboard**

88 character keyboard allows user to enter IDs directly into CFEZ.



A13-1014-254

**USB Barcode Reader**

Barcode Scanner scans product IDs directly into CFEZ. Automatically detected at USB port.



A13-1014-259

**USB Printer**

USB printer allows hard copy measurements to be printed. requires a standard USB cable, sold separately. Uses Thermal Printer paper, sold separately.

A21-1013-859

**USB Cable**

USB cable for use with CFEZ and USB Printer and CFEZ.

L02-1014-440

**Thermal Printer Paper**

Thermal printer paper for use with A13-1014-259 USB printer. 5 rolls per package.

For more information, please contact HunterLab at 703-471-6870, [sales@hunterlab.com](mailto:sales@hunterlab.com) or visit [www.hunterlab.com](http://www.hunterlab.com)