

# The World's Most Trusted Name In Color Quality

# HunterLab

## Introduces

## The ColorFlex EZ




### 45/0 Design: For Relentless Perfection in Color Quality

### The Power to See Color the Way Your Customers Do

HunterLab's next generation Color Flex EZ spectrophotometer takes color quality control to its highest level with 45/0 design for the ultimate in color measurement preciseness. Combining versatility, simplicity and performance, the ColorFlex EZ reflects almost 60 years of color measurement innovation in one easy-to-use, compact instrument from the world's most trusted experts in color quality.

HunterLab's ColorFlex EZ offers the advantage of 45/0 design ensuring that you not only get your colors right - but perfect every time. It is the optimal and only instrument design that measures your samples the way the human eye does. It lets you see your colors exactly the way your customers do, not just in the lab but in the real world.



**HunterLab**  
Measure Color...Measure Quality

## More Connectivity


More options to connect to a wide variety of devices:

- USB flash drives to transfer results to several machines or computers
- Keyboards to enter custom sample Id's manually
- Barcode scanners to scan samples fast without data entry errors
- Printers to get instant hard copies of sample readouts
- Computers to plot data using HunterLab's Easy Match QC software

## More Simplicity

Easy to operate for all your color measurement challenges:

- Ergonomic push buttons to deliver lightning fast results
- A large 3.7" [93.98mm] high resolution easy-on-the-eye color display
- 3 USB 2.0 communications ports to connect to a wide variety of devices
- Small size requiring very little benchtop space
- Versatility to inspect raw materials or finished products in the lab or in production
- A variety of sample handling accessories



## The HunterLab Advantage

The ColorFlex EZ is backed by over 5 decades of color quality innovation and experience from HunterLab, the world's most trusted color quality experts.

With an unmatched reputation for delivering the right solution for the right challenge, HunterLab tailors products and technologies for every color measurement need and budget offering the broadest range of color measurement solutions in the industry.

## More Performance

More capabilities for increased performance and better outcomes:

- Sampling versatility to measure everything from opaque solids, liquids, powders, granules and pellets to translucent solids and liquids
- Easy data transfer to move set ups and standards from one instrument to another and up load or download to USB Flash drives for back ups
- Expanded sample storage for 250 setups, product standards and tolerances and up to 2000 sample measurements



## More Flexibility

Built-in software using 250 set up configurations with product standard values and tolerances for a wide variety of parameters including:

- Color Scale, Indices, Illuminant, Observer to specify conditions under which data is reported
- Standard values to identify product standards or targets for color comparisons
- Pass fail tolerances to measure against product specifications
- Averaging to provide accurate measurements of uneven or inconsistent samples
- Display formats to report measurement data multiple ways such as color data, color difference data, spectral difference data and many other ways.

