

# Accessories



## ColorFlex®



- 04-6623-10** 1.25 inch (32mm) (45°/0° only)  
**02-6624-05** 1.25 inch (32mm) with glass (45°/0° only)  
**D02-1010-618** 1.25 inch (32mm) with 420nm UV Filter Assembly (45°/0°)  
**04-6623-09** 1 inch (25mm) (45°/0° only)  
**02-6624-01** 1 inch (25mm) with glass (45°/0° only)  
**04-6623-08** 0.75 inch (19mm) (45°/0° only)  
**04-6623-05** 0.5 inch (13mm) (45°/0° only)

### Port Inserts

**04-6622-00**

### Port Insert for Sample Cup


2.5 inch (64mm, 45°/0° only) The use of the sample cup enables the measurement of powders, granules, pellets, pastes and liquids. The port insert is designed to hold the sample cup in proper position at the measuring port.



**04-7209-00**

### Glass Sample Cup

2.5 inch (64mm)



**02-4579-00**

### Ring and Disk Set

For use in the glass sample cup, used for transparent and translucent liquids on semi-solids where path length must be fixed. Does not include sample cup.



**04-4000-00**

### Sample Cup Opaque Cover

Provides a light trap to exclude the interference of external light on the sample cup.



D02-1010-459

**Sample Clamp / Port Forward Stand**

Securely holds the instrument in a port forward orientation.  
The sample clamp permits to be held flat against the sample port.



D02-1010-619

**Hard-Sided Carrying Case**



D02-1010-327

**Footswitch Assembly**

The footswitch is used to initiate sample measurement without the use of hands.



L01-1010-678

**External Printer Kit, 110V**

Includes one spare ribbon, one spare roll paper and cable

L01-1010-679

**External Printer Kit, 220V**

Includes one spare ribbon, one spare roll paper and cable

A13-1010-676

**Ribbon, Printer (spare)**

D02-1011-833

**Paper, Printer, 2.25 inch wide (spare)**

For more information, please contact HunterLab at 703-471-6870, [sales@hunterlab.com](mailto:sales@hunterlab.com) or visit [www.hunterlab.com](http://www.hunterlab.com)