

Accessories



D25LT Colorimeter


For D25LT with Sample Port Down



D04-1013-402

Sample Dish

Anodized aluminum dish, 7 inch (17.8cm) diameter, used with port down sample presenter.



D04-1013-411

Sample Dish with Mount for a Petri Dish

Anodized aluminum dish, 7 inch (17.8cm) diameter, with positioning holder to support a Petri dish. Used for measuring powders, granules and flakes, such as spices.


For D25LT with Sample Port Up



04-7767-00

Glass Sample Cup

A 5 inch (12.7cm) sample cup that is ideal for the analysis of large granular products and foods.



12-0001-20

Cover Glass

A 4 inch (10.2cm) square cover glass to prevent dust and sample material from contaminating the inside of the sensor.



L02-1013-431

Sample Dish Assembly and Polycarbonate Sensor Window

Provides an anodized aluminum sample dish with polycarbonate bottom, a lid for the sample dish, a polycarbonate sensor port window (not shown), and a black light trap to be used in place of the black glass instrument standard.

L02-1008-437

Sample Dish with Polycarbonate Bottom

Round, 7 inch (17.8cm) diameter, anodized aluminum sample dish with polycarbonate bottom. One included with L02-1013-431.

A04-1001-308 Polycarbonate bottom (replacement)

L02-1008-436

Sample Dish with Glass Bottom

Round, 7 inch (17.8cm) diameter, anodized aluminum sample dish with glass bottom. Glass bottom is more scratch resistant than polycarbonate, but can not be used in food processing areas.

A04-1002-920 Glass bottom (replacement)

L02-1008-435

Lid for Sample Dish

Round lid, with handle, used to cover sample dish. One included with L02-1013-431.

A02-1002-965

Light Trap

A round, black light trap, that is to be used in place of the black glass instrument standard when using sample dishes. One included with L02-1013-431.

D02-1013-291

Measurement Port Polycarbonate Window

A polycarbonate window that can be mounted in the D25LT measurement port. Prevents dust and sample material from contaminating the inside of the sensor. One included with L02-1013-431.

For All D25LT

L02-1013-427

D25LT External Printer

External printer with USB input.

A10-1013-423

USB 2.0 Flash Drive, 256MB

Enables software updates, data transfer, and driver installations.

A10-1013-424

USB 2.0 Port Hub

Provides 4 USB ports for connection with other USB supported hardware.

A10-1011-777

USB-To-Serial Adapter

For more information, please contact HunterLab at 703-471-6870, sales@hunterlab.com or visit www.hunterlab.com