

ColorFlex® EZ Citrus

measure COIOT the way your COIOT EVE sees it



The ColorFlex® EZ Citrus

Like all members of the ColorFlex® EZ family of products, the ColorFlex® EZ Citrus is specifically built to meet and exceed your color measurement needs. Every ColorFlex® EZ Citrus comes with 60 years of color quality control experience and innovation by HunterLab,

The world's true measure of color.

Using a **45°/0° geometry**, the easy-to-use ColorFlex® EZ Citrus provides you with the highest level of color quality control. More importantly, it allows you to

measure COIOT the way your **EYE sees it**

The ColorFlex® EZ Citrus is specifically designed to measure the color of fresh and concentrated orange, grapefruit, and lemon juice.

Citrus color scales provided in firmware include:

- Citrus Number
- Citrus Redness
- Citrus Yellowness



Other important color scales, including **Hunter L,a,b** and **CIE L*a*b*** are provided within the firmware, allowing for the capability to use your ColorFlex® EZ Citrus to measure the reflected color of a wide range of non-citrus samples - liquids, semi-solids, powders, and solids.



As with all ColorFlex® EZs,

the ColorFlex® EZ Citrus

offers you more...

- Simplicity with one-touch sample measurement
- Connectivity with 3 USB 2.0 ports
- Ruggedness with a sealed optical design

for the industrial environment

- Versatility standard with over 250 product set-ups, including a variety of parameters such as:
- Color Scale, Indices, Illuminant and Observer measurement condition settings
- Standards values
- Pass /Fail tolerances
- Averaging for more repeatable measurement of non-uniform samples



Inside its optional, custom designed, watertight carrying case will be everything you will need to immediately start measuring citrus color values...

- ColorFlex[®] EZ Citrus instrument
- Complete set of instruments and diagnostic standards
- Calibrated Orange Juice 4 Standard
- Orange Juice Tube Holder
- Sample Clamp / Port Forward Stand

... all backed up by the finest customer support in the industry.



For more information go to www.hunterlab.com or contact your local HunterLab representative.





Accessories



A11-1011-510

Orange Juice Standards, Set of 6 (OJ1-OJ6)

Plastic and uncalibrated tubes with positioning screw. It is used for hitching the instrument to USDA's standards



D02-1011-127

Port Insert: Orange Juice Tube Holder

25 mm (1 inch) in diameter orange juice tube holder. It is used to hold sample tubes at the measurement port.

D02-1010-459

Sample Clamp / Port Forward Stand

Securely holds the instrument in a port forward orientation. The sample clamp permits samples to be held flat against the sample port.





A13-1014-294

A13-1014-254

A13-1014-259

USB Flexible Keyboard

88 character keyboard allows user to enter IDs directly into ColorFlex® EZ Citrus. Sold separately.

USB Barcode Reader

Barcode Scanner scans product IDs directly into ColorFlex® EZ. Automatically detected at USB port. Sold separately.

USB Printer

USB printer allows hard copy measurements to be printed. Requires a standard USB cable, sold separately. Uses Thermal Printer paper, sold separately.



ISO 9001 Certified; **CE**Certified Hunter Associates Laboratory Inc., 11491 Sunset Hills Road, Reston, VA 20190-5280 Tel 703.471.6870 • Fax 703.471.4237 • sales@hunterlab.com • www.hunterlab.com